使用说明书	中文	页码	1
FireMagic™ 3D 冷盐水灌注射频消融导管			
Instructions for Use	English	Page	7
FireMagic™ 3D Irrigated Ablation Catheter			
Kullanım Talimatları	Türkçe	Sayfa	12
FireMagic™ 3D Sulu Ablasyon Kateteri			
Οδηγίες χρήσης	Ελληνικά	Σελίδα	17
FireMagic™ Καθετήρας κατάλυσης 3D με λειτουργία έκπλυσης			
Instrucciones de uso	Español	Página	23
Catéter de ablación irrigado 3D FireMagic™			
Instruções de Utilização	Português	Página	29
Cateter de Ablação com Irrigação 3D FireMagic™			
Инструкция за употреба	Български	Страница	35
FireMagic™ 3D иригиращ абляционен катетър			
Инструкция по применению	Русский	Страница	41
Катетер абляционный 3D орошаемый FireMagic™			
Bedienungshinweise	Deutsch	Seite	47
FireMagic™ 3D gekühlter Ablationskatheter			
Instructions d'utilisation			53
Cathéter d'ablation irrigué 3D FireMagic™	Français	Page	
Istruzioni per l'uso			59
FireMagic™ Catetere Ablatore 3D Irrigato	Italiano	Pagina	
Instructies voor het gebruik			65
FireMagic™ 3D geïrrigeerde ablatiekatheter	Nederlands	Pagina	
Інструкція з застосування	Українська	Сторінка	71
FireMagic™ 3D Абляційний катетер зі зрошенням			

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FireMagic™ 3D Irrigated Ablation Catheter

Shanghai MicroPort EP MedTech Co., Ltd.

【CATHETER DESCRIPTION】


Fig. 1: FireMagic™ 3D Irrigated Ablation Catheter

The FireMagic™ 3D Irrigated Ablation Catheter (FireMagic 3D for short) has a high-torque shaft with a deflectable tip section containing an array of platinum electrodes. Tip deflection is controlled at the proximal end by a handle.

All electrodes can be used for the purposes of cardiac electrophysiological mapping and stimulation, and the tip electrode can also be used for ablation. For ablation, the catheter is used in conjunction with a radiorequency (RF) ablation device and dispersive pad (indifferent electrode), and the tip electrode serves to deliver RF current from the RF generator to the desired ablation site. The temperature sensor embedded in the tip electrode is used for monitoring the tip temperature. At the proximal end of the catheter, a saline input port with a standard luer servers to permit the injection of normal saline to irrigate the tip eletrode. With connection to the infusion pump, normal saline passes through the port to irrigate and cool the ablation site.

When used with a 3D cardiac electrophysiological mapping system, FireMagic 3D can detect the location of the tip electrode via a magnetic location sensor embedded in the catheter, and can conduct 3D cardiac modeling in Columbus™ 3D EP Navigation System. The ECG and ablation point are traced on the established 3D cardiac model to guide the analysis of abnormal location of arrhythmia and provide accurate ablation. By embedding another magnetic location sensor in the end of the deflectable tip section, Columbus™ 3D EP Navigation System allows real-time display of the catheter curvature.

The FireMagic 3D is an active invasive device for patients aged 18 or older. The duration of use is short term. Guided by a transseptal sheath, the device is pushed along the vein to the target site of the cardiac chamber. The catheter is designed for easy tip manipulation in order to reach the desired ablation sites. Throuth the tip ablation electrode, radiofrequency power may be delivered to the lesion sites to achieve the therapeutic purpose. According to the requirement in IEC60601-1, the catheter and its connection cable are both defined as applied part.

【PRODUCT SPECIFICATION】

The FireMagic 3D has a color code near the handpiece. Different code colors show different shapes of the curve. For specifications and dimensions of the product available, see Table 1. Doctors can make choice in accordance with the location of disease and their own personal habits.

Table 1: Specification for FireMagic™ 3D Irrigated Ablation Catheter

Types	Shape of curve /Color	Length (cm)	Size* (F)	Electrode interval (mm)	Length of tip electrode (mm)
EPN8BTC	B/red	105	8	2-5-2	4
EPN8DTC	D/blue	105	8	2-5-2	4
EPN8FTC	F/orange	105	8	2-5-2	4
EPN8JTC	J/black	105	8	2-5-2	4

*8F=2.67mm

【INTENDED USE】

The FireMagic™ 3D Irrigated Ablation Catheter is indicated for catheter-based cardiac electrophysiological mapping (stimulating and recording) in conjunction with a 3D EP navigation system, and when used with the RF generator, for the treatment of atrial fibrillation (AF) via catheter-based radiofrequency ablation in patients of 18 or older age.

【INDICATIONS】

Drug refractory symptomatic paroxysmal atrial fibrillation and persistent atrial fibrillation

【CONTRAINDICATIONS】

- Left atrial anterior and posterior diameter $\geq 50\text{mm}$ or left ventricular ejection fraction $<40\%$;
- Patients with advanced heart failure (NYHA Functional Class III–IV);
- Cerebral vascular diseases including apoplexy and transient ischemic attack within the last two weeks;
- Acute myocardial infarction within 3 months;
- Unstable angina;
- Patients with intracardiac mural thrombus or subjected to ventriculotomy or atriotomy in the past 4 weeks;
- The catheter may not be appropriate for patients with artificial valve, whereas it can be used for patients with aortic heart valve if no retrograde approach through aortic heart valve is performed;
- General active infection;
- For patients with left atrium thrombus (LATH), mucous tumor or interatrial septum or patch, septum perforation is not advised;
- For patients who have been subjected to aortic valve replacement, retrograde transaortic access procedure is not advised;
- Patients with obvious bleeding liability and hematologic disorder;
- Use in children.

【TARGET GROUP】

For 18-75 years old adults, men or non-pregnant women

【WARNING】

- With significant potential X-ray exposure, catheter ablation may subject patients and operators to acute radiation damage and increase the risk of

impact on their bodies and genes. Catheter ablation can be applied only when sufficient attention is paid to the potential exposure involved in the procedure and measures are taken to minimize radiation exposure.

- Implantable pacemakers and implantable cardioverter/defibrillator (ICDs) may be adversely affected by RF current. It is important to: have temporary external sources of pacing and defibrillation available during ablation; temporarily reprogram the pacing system to minimum output or OFF mode to minimize the risk of inappropriate pacing. Exercise extreme caution during ablation when in close proximity to atrial or ventricular permanent leads; program the ICD to the OFF mode during the ablation procedure; and, perform complete implantable devices analysis on all patients after ablation.
- The temperature sensor in the catheter electrode is to measure the electrode temperature instead of the tissue temperature. Because of the cooling effect that normal saline infusion has on the electrode and the feature that RF ablation is passively heated, the electrode temperature displayed on RF generator is noticeably lower than that of the tissue in the ablation zone during the procedure. Temperature sensor is used to check if the infusion flow is appropriate, and at the same ablation temperature, larger infusion flow corresponds to higher RF discharge power and generates a deeper ablation zone. Therefore, special attention shall be paid to the infusion flow and actual power of ablation in procedure.
- Throughout the procedure, balance of patient's body fluid needs to be monitored to prevent overload of body fluid. Some patients, whose capability to handle overload of body fluid drops for some reasons, tend to suffer from acute pulmonary edema or heart failure during or after the procedure. This is particularly true of patients with congestive heart failure and elderly people. It is important to evaluate the level of risk that patients may develop overload of body fluid.
- Performing procedure again with this catheter on patients who has been subjected to RF ablation may face higher risk of perforation or pericardial effusion.
- When using 3D Irrigated ablation catheter, ferromagnetic devices and equipment will have potential offset interference with the product in picking up the signal for displaying ablation electrode and curve position. Therefore, throughout the procedure from point acquisition for 3D reconstruction to completion of ablation, approaching to or moving ferromagnetic devices and equipment such as MRI equipment and X-ray apparatus during procedure shall be avoided. In case an offset is present, usually returning the device to the original location will remove the offset.
- So far, no proof shows that this device remains safe at $>60^{\circ}\text{C}$.
- Never immerse the near-end handle or cable connector in liquid. Otherwise, the electric performance may be affected.
- This product shall be operated by trained personnel in electrophysiologically well-equipped catheter room.
- In insertion, never exert excessive force, and pull out the catheter in case of resistance.
- In cleaning, never wipe or twist the ablation end electrode heavily.
- To ensure patient's safety and the integrity of electrode catheter, never wipe the catheter with alcohol.
- Make sure that the output voltage used by the device is below rated accessory voltage of the catheter and connection cable.
- When placing the catheter and connection cable, make sure that they are not in contact with the patient or other cables. When they are not in service for the time being, keep them in a place isolated from the patient.
- When used, catheter shall be disposed of as medical waste according local laws and regulations.
- When the catheter is connected to the supply mains-operated stimulator, inadvertently introducing leakage current may enter into the heart. Please use the internal power -operated simulator.
- Before handling the external pulse generator, the patient cable or indwelling leads, steps shall be taken to equalize the electrostatic potential between the user and the patient, for example by touching the patient at a site remote from the pacing lead.
- The catheter shall be connected to the non-implantable pulse generator before the pacing leads are connected to the catheter.
- When handling indwelling leads, the terminal pins or exposed metal are not to be touched nor be allowed to contact electrically conductive or wet surfaces.
- This device is packaged for single use only. Do not reuse, reprocess or resterilize. Reuse, reprocessing or resterilization may compromise the structural integrity of the device and/or lead to device failure that in turn may result in patient injury, illness or death. Also, reprocessing or resterilization of single use devices may create a risk of contamination and/or cause patient infection or cross-infection, including, but not limited to, the transmission of infectious disease(s) from one patient to another. Contamination of the device may lead to injury, illness, or death of the patient.
- If the patient with a history of heparin induced thrombocytopenia (HIT), bivalirudin is a preferred alternative to heparin for anticoagulation during procedures.

【PRECAUTIONS】

- In case of obvious low power output, readings of high impedance, abnormal signal disturbing and the like, immediately check if the dispersive pad is stuck and damaged, the infusion tube unblocked, the external reference patch securely bonded, and the connecting cable reliably connected; and if the abnormality cannot be removed, do not increase the ablation power.
- To avoid thromboembolism, usually intravenous heparin infusion should be given to maintain an activated clotting time >300 s during ablation or anticoagulation medicine is taken for about 3 months after the procedure as instructed by the doctor, giving NOAC priority. No consensus is available yet as to the necessity for short-time anticoagulation after ablation.
- Flush the catheter and infusion tube. Before pushing the catheter into the body, perform complete degassing to prevent air embolism, and during the service, always maintain certain low infusion flow rate to prevent the blood from coagulation in the infusion tube.
- Make sure to carefully follow the power-adjusting steps specified in the

operating instructions. Excessively increasing the ablation power may cause steam burst and perforation.

- While operating the catheter, always be careful to avoid damage to the heart which causes perforation or tamponade.
- Before pushing in or pulling out the catheter, always fully pull back the thumbknob so as to straighten the catheter.
- It is recommended that the ablation operations for each catheter be less than 50 times, with 2 minutes for each operation.
- FireMagic 3D used in conjunction with RF ablation device can transfer considerable electric energy. Improperly operating the catheter and dispersing electrode (especially when operating the device) may damage the patient or operator. During the energy transfer, patient must not touch any earthed metal surface. If the tip temperature does not rise during ablation, immediately stop ablation and replace the catheter.
- Regularly check and test the reusable Connection cable and accessories.
- For connection, use equipment and cables specified in the operating instructions and compatible with this device.
- Prior to use, carefully read the operating instructions for the product. Pay attention to all warnings and precautions in the instructions. Otherwise, accident may occur.
- Do not use near MRI equipment since movement or heating of the catheter may occur and the image on the display may become distorted.

【ADVERSE REACTIONS】

Pulmonary embolism, myocardial infarction, cardiac perforation/tamponade, vascular bleeding/local hematomas, thrombosis, AV fistula, pseudoaneurysm, thromboembolism and vasovagal reactions, pericardial effusion/tamponade, thrombi, air embolism, arrhythmias and valvular damage, pneumothorax and hemothorax, pulmonary edema, hypoxia, pleural effusion, acute respiratory distress syndrome(ARDS), congestive heart failure, aspiration pneumonia, pneumonia, asthmatic attack, hypotension, pacemaker/ICD malfunction, anemia, thrombocytopenia, disseminated intravascular coagulation, epistaxis, systemic infection, urinary tract infection, sedation induced apnea, sedation induced CO2 retention with lethargy and cholecystitis, chest pain/discomfort, ventricular tachyarrhythmia, stroke/TIA, CVA, complete heart block, coronary artery spasm, coronary artery thrombosis, coronary artery dissection, cardiac thromboembolism, pericarditis, valvular damage and increased phosphokinase level, urinary retention, transient extremity numbness, parkinson's disease and gastrointestinal diverticulosis, pulmonary vein stenosis, atrio-oesophageal fistula, silent strokes, death, electrical shock, open-heart surgery, prolonged procedure, scratches of intravascular wall or myocardial wall

【DEVICES USED IN COMBINATION】

- Transeptal Sheath matching the catheter, 8.5F
- Radiofrequency (RF) Generator
- Cable for the 3D Irrigated Ablation Catheter
- Columbus™ 3D EP Navigation System
- Irrigation Pump

【CONNECTION FOR PROCEDURE】


Fig. 2: The sketch for connection

4. FireMagic 3D 5. Catheter handle 6. Strap for irrigation pump 7. Standard Luer
 8. Cable for the 3D Irrigated ablation catheter 9. Connection cable
 10. Cable for the 3D irrigated ablation catheter 11. Columbus™ external reference patch
 12. Connector for the external reference patch

As shown in Fig. 2, connect FireMagic™ 3D Irrigated Ablation Catheter with relevant procedure equipment. Equipment compatible or used in combination include:

Compatible Equipment	Brand and Type	Remarks
RF Generator	Stockert EP Shuttle	It should be connected through Columbus™ 3D EP Navigation System
	HL-100F	
	IBI 1500T	
	Atakr II	
Irrigation Pump	Coolflow	/
3D EP Navigation System	Columbus™	/
External Reference Patch		

Note: 1) The MAX output frequency for compatible RF generator should be at the range of 450KHZ-575KHZ

2) Please choose the neutral electrode which is recommended by the manufacturers of the compatible RF Generator, for the details, refer to the IFU of the generator, if not, please choose the neutral electrode obtained the CE approve.

• Connection Cable

The cable is chosen and used for connection of FireMagic 3D with 3D electrophysiological mapping system. Via connection cable, the FireMagic 3D can

interface with external equipment. The connection cable is the accessory for FireMagic 3D catheter.

Name	Types
Cable for 3D Irrigated ablation catheter	EPS 150 C

Note: The cable is the product of MicroPort EP.

The rated accessory voltage of catheter and connection cable is 450V

• **Connection with Columbus™ 3D EP Navigation System**

Select specified type of cable. Based on the corresponding “key” position, plug the plug on the side marked “catheter” of the connection cable in the jack in the catheter, and a “click” indicates the completion of the plugging; and plug the plug on the other side of the cable in the jack marked “Map interface” in the PIU of Columbus™ 3D EP Navigation System based also on the corresponding “key” position and a “click” indicates the completion of the plugging.

Remove the cable from the connecting devices on the two sides: just pull it out with slight force by holding the arrow on the cable connector.

• **Connection with irrigation pump**

Butt and tighten the luer on the compatible irrigation pump strap to the saline tube luer on the handle of FireMagic 3D.

• **Connection with other devices and equipment**

1) For the connection and use of external reference patch, refer to IFU for Columbus™ External Reference Patch;

2) For the connection with and use of Columbus™ 3D EP Navigation System, refer to IFU for Columbus™ 3D EP Navigation System.

【SUGGESTED INSTRUCTIONS FOR USE】

- To choose a catheter, note that the shape of the curve suits the heart anatomy to be reached;
- First, carefully check if the package of the catheter is intact;
- Remove the catheter from the package, carefully check the integrity of the electrodes and catheter, and place them on the sterile bench;
- Establish channel in the peripheral great vessels by sterile technology;
- Screw the catheter luer to the strap connector of the infusion pump, vent by opening to the max infusion flow rate ($> 30\text{ml/min}$) and check that the infusion channel is unblocked: all the 6 saline irrigation holes in the ablation electrode bleed water normally;
- Switch the Irrigation pump flow rate to and maintain it at the low flow rate without ablation (2-5 ml/min), and insert the catheter into the established vessel channel;
- Connect the catheter to the RF generator or 3D equipment with compatible connecting cable, and check that the connection with the RF generator, multi-channel electrophysiological recorder, stimulator and the like is properly performed;
- Push the catheter into the endocardium zone to be examined and use X-ray and ECG to help with correct locating;
- Change the bending of the tip with the thumbknob to help locating by the catheter tip. When the thumbknob is pushed forward, the tip is deflected. When the thumbknob is pulled back, the tip straightens;
- Operate the catheter so that the tip electrode arrives at the target location in the heart cavity, perform point acquisition for 3D reconstruction (when using the FireMagic™ and connecting the Columbus™ system) or pick up ECG;
When the position of the tip electrode or catheter bending cannot be displayed or the data is inaccurate, check if the connecting interfaces are correctly stabilized; or replace the cable or catheter. If the replacement results in returning to normal (the normal condition), it indicates the damage to the device;
- After selecting the target location for ablation, set the irrigation pump flow rate and ablation parameters, switch to high irrigation flow rate and start ablation; Recommendation for selection of infusion flow: prior to ablation, infuse at high flow for 5s, and maintain the high flow within 5s at the completion of ablation; for an ablation power of 15-30W, use a high flow at 17ml/min, and for an ablation power of 30-50W, use a high flow at 30ml/min.
Before applying RF energy, verify that the infusion flow is increased by lowering the ablation electrode temperature;
And the duration for each ablation shall be within 120s. During ablation, it is allowed to move the ablation electrode to the next location;
When starting the RF current, the circuit impedance shall be about 100ohms; During discharge, carefully check the output power. An excessively low output power may suggest that the Saline irrigation holes in the ablation electrode are blocked. At that time, discharge must stop and the catheter must be pulled out. Clean the coagulum in the ablation electrode, verify that the Saline irrigation holes bleed water normally, and re-insert after complete venting is done: clean the tip gently with sterile gauze pad soaked with sterile saline; and never rub or twist the ablation electrode heavily. Otherwise, the bonding of the ablation electrode may be damaged, which leads to the electrode loosening. Alternatively, flush the infusion tube by injecting normal saline at high pressure with a syringe.
- At the completion of the procedure, remove the catheter and dispose of it as specified.
- In case of any questions on usage or product performance, consult local distributor or the manufacturer.

【TRANSPORTATION REQUIREMENTS】

In transit, the product shall be protected from heavy load, direct sunlight and rain or as specified in the ordering contract. The temperature during transportation shall be kept between 0℃ and 45℃.

【STORAGE REQUIREMENTS】

The product shall be stored in a shady and cool, dry, clean, and well-ventilated warehouse which is in natural air circulation environment. The temperature during storage shall be kept between 0℃ and 45℃.

【SHELF LIFE】

The shelf life of product is three years while it meets the conditions for storage.

【STERILIZATION】

This product has been sterilized with ethylene oxide gas. Never re-sterilize and reuse it. Do not use the catheter if the package is open or damaged. Use the catheter prior to the expiration date shown on the package label.

【SYMBOLS EXPLANATION】

1.		DO NOT REUSE
2.		CONSULT INSTRUCTIONS FOR USE
3.		PROTECT FROM HEAT SOURCE AND RADIATION SOURCE
4.		KEEP DRY
5.		USE BY
6.		BATCH CODE
7.		CATALOGUE NUMBER
8.		STERILIZED USING ETHYLENE OXIDE
9.		DO NOT USE IF PACKAGE IS DAMAGED
10.		QUANTITY OF PRODUCT CONTAINED 1
11.		TEMPERATURE LIMITATION
12.		TYPE OF APPLIED PART
13.		DATE OF MANUFACTURE
14.		MANUFACTURER
15.		AUTHORISED REPRESENTATIVE IN THE EUROPEAN COMMUNITY
16.		IMPORTER

【AFTER-SALES SERVICE】

With “providing the medical sector with high quality and efficacy medical products” as its top operational objective, Shanghai MicroPort EP MedTech Co., Ltd. (hereinafter referred to as MicroPort EP Co.) guarantees that its products are free of defects in materials or manufacturing when the clients receive them. For other questions relating to the products, directly consult the company.

【SOLEMN STATEMENT】

Shanghai MicroPort EP MedTech Co., Ltd. expressly states herein that its **FireMagic™ 3D Irrigated Ablation Catheter** is disposable and cannot be reused. MicroPort EP Co. will not recommend, indicate and imply in any manner the reusability of the system, and will not assume the responsibility for any accident or product damage resulting from reuse.

FireMagic™ 3D Irrigated Ablation Catheter can be connected and used only with the compatible devices specified herein, and MicroPort EP Co. will not assume the responsibility for damage to product device, procedure failure and the like resulting from operating mistakes or any other human factors.

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Shanghai MicroPort EP MedTech Co., Ltd.

Address: Building 23&26&28, Lane 588, Tianxiong Rd., 201318 Shanghai, PEOPLE'S REPUBLIC OF CHINA

Postcode: 201318

Tel: +86 (21) 60969600

Fax: +86 (21) 20903925

E-mail: customerservice@everpace.com

Website: www.everpace.com

【AUTHORIZED REPRESENTATIVE IN THE EUROPEAN COMMUNITY】

MicroPort Medical B.V.

Address: Paasheuvelweg 25, 1105BP Amsterdam, The Netherlands

E-mail: cs@microport-int.com

Tel: +31(0)205450100-8

Fax: +31(0)205450109


The copyright of this document is owned by Shanghai MicroPort EP MedTech Co., Ltd.

Copying is not allowed without its consent.

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**De copyright van dit document is eigendom van Shanghai MicroPort EP MedTech