OptimAblate™

使用说明书	中文	页码	1
灌注泵管系统			
Instructions for Use	English	Page	4
OptimAblate™ Tubing Set			
Kullanım Talimatları	Türkçe	Sayfa	7
OptimAblate™ Tüp Seti			
Οδηγίες χρήσης	Ελληνικά	Σελίδα	10
OptimAblate™ Σετ σωλήνων			
Instrucciones de uso	Español	Página	13
Set de tubos OptimAblate™			
Instruções de Utilização	Português	Página	16
Conjunto de Tubos OptimAblate™			
Инструкция за употреба	Български	Страница	19
OptimAblate™ система за инфузии			
Инструкция по применению	Русский	Страница	22
Комплект трубок OptimAblate™			
Gebrauchsanweisung	Deutsch	Seite	25
OptimAblate™ Schlauchset			
Instructions d'utilisation	Français	Page	28
Ensemble de tubes OptimAblate™			
Istruzioni per l'uso	Italiano	Pagina	31
OptimAblate™ set di tubi			
Instructies voor het gebruik	Nederlands	Pagina	34
OptimAblate™ slangenset			
Інструкція з застосування	Українська	Сторінка	37
OptimAblate™ Набор трубок			

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Shanghai MicroPort EP MedTech Co., Ltd.

【WARNING】

- Do not expose the tubing set or connecting pieces to organic solvents such as isopropyl alcohol (IPA).
- Prior to use, carefully check the product shelf life on the package and do not use product with damaged package and beyond the shelf life.
- This device is intended for one time use only. Do not reprocess or reuse. Reuse can cause patient injury and/or the communication of infectious disease(s).
- This tubing set is intended for use with the OptimAblate™ Irrigation Pump only. Refer to the pump Operator Manual for additional Warnings and Precautions.
- When used, this device shall be disposed of as medical waste according to local laws and regulations.

【PACKING LIST】

- A OptimAblate™ Tubing Set
- An instructions for use

【OVERVIEW】

The OptimAblate™ Tubing Set is a special single use accessory of the OptimAblate™ Irrigation Pump. The OptimAblate™ Tubing Set cooperate with the OptimAblate™ Irrigation Pump, linkage with cardiac RF generator, both ends connect with the infusion bag (or infusion bottle) and the irrigated ablation catheter, and the normal saline at room temperature can be accurately irrigated at constant flow rate to the irrigated ablation catheter tip, which plays a role in local tissue cooling and preventing thrombosis, improving the ablation efficiency and achieving to improve the operation success rate.

【PRODUCT SPECIFICATION】

Table 1: OptimAblate™ Tubing Set Specification

No.	Model	Effective length
1	EPAB 3200	3200 mm

【INDICATIONS FOR USE】

OptimAblate™ Tubing Set is cooperatively used with OptimAblate™ Irrigation Pump only, which provides irrigation solutions.

【CONTRAINDICATIONS】

None.

【ACCESSORIES REQUIRED FOR USING PRODUCT】

Quantity	Product	Model	Manufacturer
1	OptimAblate™ Irrigation Pump	EPE-IGP-1A	Shanghai MicroPort EP MedTech Co., Ltd.

【RISKS AND SIDE EFFECTS】

No found risk and side effects.

【SUGGESTED INSTRUCTIONS FOR USE】

- 1) Carefully inspect the outside of the box for indications of damage, and then remove the sealed package from box. Do not use if the package is open or damaged.
- 2) Inspect the tubing set carefully for integrity and overall condition. Do not use if the tubing is kinked.
- 3) Remove the tubing set from its package using standard safe hospital practices (aseptic technique).
- 4) Make sure the 3-way stopcock is in the closed position prior to use. Ensure that all connecting pieces in the tubing set are secure.
- 5) Ensure the irrigation fluid is at room temperature. Fully spike the container and hang it on the pole above the pump. Gently squeeze and release the drip chamber until approximately ½ full.
- NOTE: Vent should remain in closed position unless solution is in a solid stopper glass bottle.
- 6) Purge air from the system by opening the 3-way stopcock and filling the tubing set with the irrigation fluid. Remove any trapped air and then close the stopcock.
- 7) Dry the outer surfaces of the tubing and bubble detectors.

NOTE: To ensure proper operation of the bubble detectors, the outer surface of the tubing must be dry.

- 8) Open the pump door and press the pump tubing against the back wall of the first bubble detector. Lock the guide key into position. Note that the pump operates in a counter-clockwise direction.
- 9) Wind the tubing around the pump head and press it into the 2nd bubble detector, locking the guide key into position. Close the pump door.
- 10) Verify tubing set integrity.
- 11) Securely connect the stopcock to a catheter (do not swab the connecting ends with IPA). To continue, refer to the pump and the irrigated ablation catheter for recommended procedures and parameters.

【STERILIZATION】

OptimAblate™ Tubing Set has been sterilized by epoxy ethane before delivery.

【STORAGE & TRANSPORTATION REQUIREMENTS】

The product shall be stored in a shady and cool, dry, clean, and well-ventilated environment free of corrosive gases. The temperature during storage shall be kept in between 0℃ and 45℃.

In transit, the product shall be protected from heavy load, direct sunlight and rain or as specified in the ordering contract. The temperature during transportation shall be below 45℃.

【MANUFACTURING DATA】

See label statement for service life.

【PERIOD OF VALIDITY】

The sterilization period of validity is two years under specified storage conditions.

【SYMBOLS DESCRIPTION】

- | | | |
|-----|---|---|
| 1. |  | QUANTITY OF PRODUCT CONTAINED 1 |
| 2. |  | STERILIZED USING ETHYLENE OXIDE |
| 3. |  | CONSULT INSTRUCTIONS FOR USE |
| 4. |  | KEEP DRY |
| 5. |  | PROTECT FROM HEAT SOURCE AND RADIATION SOURCE |
| 6. |  | NO NOT USE IF PACKAGE IS DAMAGED |
| 7. |  | TEMPERATURE LIMITATION |
| 8. |  | DO NOT REUSE |
| 9. |  | NON-PYROGENIC |
| 10. |  | USE BY |
| 11. |  | DATE OF MANUFACTURE |
| 12. |  | BATCH CODE |
| 13. |  | CATALOGUE NUMBER |
| 14. |  | MANUFACTURER |
| 15. |  | AUTHORISED REPRESENTATIVE IN THE EUROPEAN COMMUNITY |
| 16. |  | IMPORTER |

【AFTER-SALES SERVICE】

With “providing the medical sector with high quality and efficacy medical products” as its top operational objective, Shanghai MicroPort EP MedTech Co., Ltd. (hereinafter referred to as MicroPort EP) guarantees that its products are free of defects in materials or manufacturing when the clients receive them. For other questions relating to the products, please directly consult the company.

【SOLEMN STATEMENT】

The potential adverse event listed in the user manual, collected by Shanghai MicroPort EP MedTech Co., Ltd. (MicroPort EP), is based on the information, data from clinical trials and markets of China and other countries or regions instead of real condition occurred. The potential adverse event is taken as the product features sufficiently warned and explained instead of the defects or potential

risks. Subject to the current scientific level and knowledge of MicroPort EP, the product may have some undetected or unknown potential risks and that cannot be considered defects or faults of the product.

MicroPort EP will not responsible for any medical cost or direct and indirect loss caused by the improper usage (including reuse), improper model choice, maloperation or other accident caused by human, illegal purchase, risks cannot be detected by current technology or potential adverse event or risks listed in this user manual, regardless of any guarantee, contract, tort or other regulations the claim is based on.

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Shanghai MicroPort EP MedTech Co., Ltd.

Address: Building 23&26&28, Lane 588, Tianxiong Rd., 201318 Shanghai, PEOPLE'S REPUBLIC OF CHINA

Postcode: 201318

Tel: +86 (21) 60969600

Fax: +86 (21) 20903925

E-mail: customerservice@everpace.com

Website: www.everpace.com

【AUTHORIZED REPRESENTATIVE IN THE EUROPEAN COMMUNITY】

MicroPort Medical B.V.

Address: Paasheuvelweg 25, 1105BP Amsterdam, The Netherlands

E-mail: cs@microport-int.com

Tel: +31 (0)20 545 0100-8

Fax: +31(0)205450109

**The copyright of this document is owned by Shanghai MicroPort EP MedTech Co., Ltd.
Copying is not allowed without its consent.**

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De copyright van dit document is eigendom van Shanghai MicroPort EP MedTech Co., Ltd.