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FORLNK™ Cable for Ablation Catheter

Shanghai MicroPort EP MedTech Co., Ltd.

【CABLE DESCRIPTION】

The cable is designed to connect a Cardiac Radiofrequency (RF) Ablation Catheter which is produced by Shanghai MicroPort EP MedTech Co., Ltd. (hereinafter called "Microport EP" for short) to the appropriate RF generator. The cable may also be reused subjected to the reuse instructions provided below. Different types of accessory cables that suitable for compatible RF generator are shown in Products Specification Table 1. Also can select and use appropriate type of accessory cable according to the relative instruction of compatible RF generator and accessory cable in the Instruction for Use Ablation Catheter of Microport EP.

The cable is suitable to connect a FireMagic™ Cardiac RF Ablation Catheter produced by Microport EP.

【PRODUCT SPECIFICATION】

Table 1: Specification of Ablation catheter--Accessory cable

Index	Cable type	Characteristic of Connector (connect to RF generator)		Length	Type of compatible RF generator
1	EPW B4	14pins		3m	IBI 1500T
2	EPW ST	10pins		3m	MicroPort EP/Stockert EP Shuttle
3	EPW AK	10pins		3m	Atakr II
4	EPWS 300SA	26pins		3m	SmartAblate
5	EPWS 300MS	9pins		3m	Meastro
6	EPWS300 WST	10pins		3m	MicroPort EP/Stockert EP Shuttle
7	EPWS300 WB4	14pins		3m	IBI 1500T /Ampere
8	EPWS300 WSA	26pins		3m	SmartAblate
9	EPWS300 WMS	9pins		3m	Meastro
10	EPWS300 WAK	10pins		3m	Atakr II

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The cable is designed to connect a Cardiac Radiofrequency (RF) Ablation Catheter to the appropriate RF generator.

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The cable should only be used by physicians and teams who are trained in electrophysiological techniques.

【WARNING】

- Use the cable prior to the "Use by" date on the package label.
- Patient or operator injury can result from improper handling of the cable.
- If the cable is used in the presence of electric equipment, noise may be induced into the cable.
- If conventional line powered equipment is used near the patient, follow electromechanical device and hospital standards to prevent extraneous electrical currents from reaching the electrophysiological equipment, catheter and heart. Such currents cause lethal arrhythmias.
- The cable should be visually inspect for any damage (i.e. kinks, nicks, crushed or elongated sections) prior to use or reuse. Do not attempt to repair any damages. In case of doubt, discard the cable and do not reuse.
- When used, cable shall be disposed of as medical waste according local laws and regulations.

【OPERATING INSTRUCTIONS】

Connection

- 1) The connector near the cable label is for RF generator, the other connector is for the catheter. Look at the relationship in the table 1, select the compatible type of the accessory cable for the RF generator.

- 2) Connect the cable plug connector with the arrow mark to the catheter handle connector into the corresponding key bit slot: Hold the catheter handle, and push the cable connectors into the catheter handle connector. When you hear a sound "keta", that means the two connectors has locked to each other.
- 3) Connect another cable plug connector to the compatible Radiofrequency as step (2).

Disconnection

- 1) Disconnect the catheter side of the cable: Holding the catheter handle; another hand holding the arrow mark of the cable connector simultaneously and you can pull it out easily. Do not pull the body of the cable.
- 2) Disconnect the RF generator side of the cable: hold the arrow mark of the connector plug of the accessory cable which is connected to the RF generator and pull it out.

【PRODUCT TROUBLE-SHOOTING】

The following points can be used as instructions for Loop system of the catheter / cable / RF generator broken off.

- RF generator temperature read error. In standby mode, most of the RF generator can display the temperature. If the temperature is too low, as error code, the temperature sensor loop may be disconnected.
- Cannot recognize the electrophysiological signals. If the recorder has a large number of electrophysiological signal interference, the reason is that the accessory cable is not connected properly.

To solve the problem of the connector and the socket pins, please deal with it in the following order:

- 1) Please verify that the settings and the type of the temperature sensor (thermistor or thermocouple) of the RF generator are matched the used catheter and accessory cable.
- 2) Please check and verify that all the connector plugs in the correct position when diagnostic/ablation/pace. Confirm the ablation connector / adapter to connect with electrophysiological monitoring well.
- 3) Pull out the plug in question, then reconnect it. Slightly pull plug and socket move, if the feedback signal is generated intermittently, there may be something wrong with the connector/adapter of the RF generator or Electrophysiological monitoring device.
- 4) If the signals are back to normal after replacing cable of Microport EP, there may be something wrong with the cable.

【GUIDE FOR STERILIZATION AND REPROCESSING】

Please refer to the table 2 as below for the reprocessing instructions.

Table 2: Reprocessing Instructions of FORLNK™ Cable for Ablation Catheter

Manufacturer: Shanghai MicroPort EP MedTech Co., Ltd.	
Device: FORLNK™ Cable for Ablation Catheter	
WARNINGS	• Do not immerse the cable connectors in fluids. • Do not reuse the cable after ten (10) uses. • The cable should be visually inspected for any damage (i.e. cuts, kinks, nicks, crushed or elongated sections) prior to use or reuse. • Do not attempt to repair any damage. In case of doubt, discard the cable and do not use. • Do not expose cable to organic solvents. • Patient or operator injury can result from improper handling or storage of the cable. To prevent damage, store the cable in a clean and secure area
Limitations on processing	The cable is designed for a maximum of 10 processing cycles.

INSTRUCTIONS	
Initial treatment at the point of use	There are no requirements for initial treatment at point of use. The cable has been sterilized by Ethylene Oxide (ETO). The sterile state is effective only when the package is intact.
Perparation before cleaning	Visually inspect the cable for any damage such as crushed or elongated sections, cuts, kinks, or nicks. If damage is observed, do not use the cable and dispose the cable as medical waste according local laws and regulations.
Cleaning: Automated	Automated cleaning of the cable is not recommended because fluid may enter the connectors and leave moisture or residue in or on the connectors. Resulting adverse effects may include improper functioning (such as noisy signals or signal degradation) or arcing of electricity between the connector pins, which can lead to patient injury

Cleaning: Manual	● For each alligator clip: Pull off the protective sleeve in direction of the junction. ● Clean and rinse with the following agents: - Cleaning agents with a pH-value 4.....10 (note manufacturer's data) - Distilled water ● Permissible cleaning method: - Manual cleaning by brushing down - Application time: 25 min - Cleaning Solution: RUHOF ENDOZIME® AW Plus W/APA - Concentration: 1:270 - Temperature: 40-45°C ● Rinse with distilled water for 3min.
Disinfection: Automated	Automated disinfection of the cable is not recommended because fluid may enter the connectors and leave moisture or residue in or on the connectors. Resulting adverse effects may include improper functioning (such as noisy signals or signal degradation) or arcing of electricity between the connector pins, which can lead to patient injury
Disinfection: Manual	● Application time: note manufacturer's data, yet do not soak more than 1 h. Disinfect with the following agent: - Mixture of 70% isopropanol and 30% water - Immersion time: 25-30 min
Drying	Dry the cable carefully at a maximum of 120°C. Drying time: 25-30 min
Maintenance, inspection and testing	Before packaging, visually inspect the cable. If visible soil is present, repeat the reprocessing instructions in this table. If there is evidence of damage, do not use the cable and dispose the cable as medical waste according to local laws and regulations.
Packaging	Coil the cable into a loop when placing it into the pouch. Do not fold the cable because doing so may break the cable. Use a pouch large enough to prevent stress to the pouch seams and to prevent excessive bending of the cable.
Sterilization	● Pack the cable in the applicable sterile container, according to the respective sterilization procedure. ● Sterilize with steam: - 5 min hold at 134°C
Storage	Follow permissible environmental conditions for storage.
Additional information	● The used product may be properly disposed of in contaminated hospital medical waste ● Never reuse the cable for more than 10 times. Reuse beyond this limit will compromise the normal function of the connection cable. MicroPort EP Co. will not assume the responsibility for the direct or subsequent adverse consequences resulting from reuse of the connection cable beyond the limit.
Manufacturer contact	For any questions relating to the cable, please directly consult the company. Please refer to the part of 【MANUFACTURER】 in this file for detailed contact information.

Note: The instructions provided above have been validated by the medical device manufacturer as being capable of preparing a medical device for reuse. It remains the responsibility of the processor to ensure that the processing as performed (using equipment, materials, and personnel in the processing facility) achieves the desired result. This requires validation and routine monitoring of the process. Likewise, any deviation by the processor from the instructions provided should be properly evaluated for effectiveness and potential adverse consequences.

【TRANSPORTATION REQUIREMENTS】

In transit, the product shall be protected from heavy load, direct sunlight and rain or as specified in the ordering contract.

【STORAGE REQUIREMENTS】

The product shall be stored in a shady and cool, dry, clean, and well-ventilated environment free of corrosive gases. The temperature during storage shall be kept in between 0°C and 45°C.

【SHELF LIFE】

While meeting the conditions for storage, the shelf life is three years.

Use this product before the end of the shelf life indicated on the package label.

【STERILIZATION】

This product has been sterilized with ethylene oxide gas. Do not use the product if the package is open or damaged. Use the product prior to the expiration date shown on the package label.

【SYMBOLS EXPLANATION】

1.

Consult Instructions for Use

2.

Protect from heat source and radiation source

3.		Keep dry
4.		Use by
5.		Batch Code
6.		Catalogue Number
7.		Sterilized using ethylene oxide
8.		Sterilized using steam or dry heat
9.		Do not use if package is damaged
10.		Quantity of product contained 1
11.		Temperature limitation
12.		Date of manufacture
13.		Manufacturer
14.		Authorized representative in the European Community
15.		Medical device
16.		Single sterile barrier system with protective packaging outside
17.		IMPORTER
18.		UNIQUE DEVICE IDENTIFIER

【AFTER-SALES SERVICE】

With “providing the medical sector with high quality and efficacy medical products” as its top operational objective, Shanghai MicroPort EP MedTech Co., Ltd. (hereinafter referred to as MicroPort EP) guarantees that its products are free of defects in materials or manufacturing when the clients receive them. For other questions relating to the products, please directly consult the company.

【SOLEMN STATEMENT】

Shanghai MicroPort EP MedTech Co., Ltd. indicates definitely that Cardiac RF ablation catheter accessory cable is used only for connecting FireMagic™ Cardiac RF ablation catheter produced by Microport EP and compatible RF generator list in the Instruction’s specification.

There is no proof to show that the accessory cable can connect to other brands of electrophysiological catheters and RF generator. Do not attempt to connect the accessory cable to other catheters or RF generator than IFU. Microport EP will not be liable for any incidental, special, or consequential loss, damage, or expense resulting, directly or indirectly, from inaccurate connecting or any other man-caused accident.

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Shanghai MicroPort EP MedTech Co., Ltd.

Address: Building 23&26&28, Lane 588, Tianxiong Rd., 201318 Shanghai, PEOPLE’S REPUBLIC OF CHINA

Postcode: 201318

Tel: +86 (21) 60969600

Fax: +86 (21) 20903925

E-mail: customerservice@everpace.com

Website: www.everpace.com

【AUTHORIZED REPRESENTATIVE IN THE EUROPEAN COMMUNITY】

MicroPort Medical B.V.

Address: Paasheuvelweg 25, 1105BP Amsterdam, The Netherlands

E-mail: cs@microport-int.com

**The copyright of this document is owned by Shanghai MicroPort EP MedTech Co., Ltd.
Copying is not allowed without its consent.**

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