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FireMagic™ Cardiac RF Ablation Catheter

Shanghai MicroPort EP MedTech Co., Ltd.

【CATHETER DESCRIPTION】

The FireMagic™ Cardiac RF Ablation Catheter (FireMagic for short) produced by Shanghai MicroPort EP MedTech Co., Ltd. is indicated for use in the diagnosis and treatment of tachyarrhythmia. The drawing below shows the catheter structure:


Fig. 1: FireMagic™ Cardiac RF Ablation Catheter

The FireMagic™ Cardiac RF Ablation Catheter is a flexible and insulative catheter with a deflectable tip comprised of thermoplastic elastomer and several platinum electrodes, and in the tip a thermocouple temperature sensor is embedded. All of the electrodes may be used for cardiac mapping and stimulation purposes; besides, the tip electrode may also be used for ablation, it's used for monopolar ablation. For ablation, this catheter is used in conjunction with a RF generator and a dispersive electrode (reference electrode), and the large tip electrode delivers radiofrequency current from the radiofrequency lesion generator to the desired ablation site. According to the requirement in IEC60601-1, the catheter and its connection cable are both defined as applied part. The temperature sensor embedded in the tip electrode is used to measure ablation temperature, and electrical connection to the sensor is made via an accessory cable.

The catheter is used for patients aged 18 or older. The catheter is rated to operate at an altitude <2000m.

【PRODUCT SPECIFICATION】

According to the curve shape and the length of deflectable section, the catheter can be divided into six different specifications, to help doctors to perform accurate positioning within different parts of the ventricle.

Depending on the hardness of the catheter shaft, the FireMagic can be divided into two major categories, the standard type and the strong support type. The strong support type may provide extra stable support ability when used to ablate near tricuspid valve area.

FireMagic has a color code near the handle. Different colors of the code show different shapes of the curve. Strong support catheter is distinguished by adding a segment of yellow sleeve.

Table 1: Specifications, model and characteristic parameters of the FireMagic

Model		Length (cm)	Diameter* (F)	Curve /color	Electrode spacing (B1-B2-B3)	Tip electrode length (mm)
Standard type	EPA 7 A TC	115	7	A/yellow	2-5-2	4
	EPA 7 B TC	115	7	B/red	2-5-2	4
	EPA 7 C TC	115	7	C/green	2-5-2	4
	EPA 7 D TC	115	7	D/blue	2-5-2	4
	EPA 7 E TC	115	7	E/white	2-5-2	4
	EPA 7 F TC	115	7	F/orange	2-5-2	4
Strong support type	EPA 8 B TC	90	7	B/red	2-5-2	4
	EPA 8 D TC	90	7	D/blue	2-5-2	4
	EPA 8 E TC	90	7	E/white	2-5-2	4
	EPA 8 F TC	90	7	F/orange	2-5-2	4

*7F=2.33mm

【INTENDED USE】

FireMagic™ Cardiac RF Ablation Catheter is indicated for catheter-based cardiac electrophysiological mapping (stimulating and recording), and when used with radiofrequency generator, for the treatment of atrioventricular nodal reentrant tachycardia and atrioventricular reentrant tachycardia via catheter-based radiofrequency ablation.

【INDICATIONS】

Atrioventricular nodal reentrant tachycardia and atrioventricular reentrant tachycardia

【CONTRAINDICATIONS】

- Patients who have intra-cardiac mural thrombus or have had a ventriculotomy or atriotomy within the preceding four weeks.

- Patients with prosthetic valves.
- Active systemic infection.
- The trans-septal approach is contraindicated in patients with left atrial thrombus or myxoma, or interatrial baffle or patch.
- The retrograde trans-aortic approach is contraindicated in patients with aortic valve replacement.
- Patients with obvious hemorrhagic tendency and have hemopathy.
- Patients with malignant tumor and have terminal disease.
- Unstable angina and acute myocardial infarction within three months.
- Stroke and transient ischemic attack within the last two weeks.

【TARGET GROUP】

For 18-75 years old adults, men or non-pregnant women

【WARNING】

- This device is packaged for single use only. Do not reuse, reprocess or resterilize. Reuse, reprocessing or resterilization may compromise the structural integrity of the device and/or lead to device failure that in turn may result in patient injury, illness or death. Also, reprocessing or resterilization of single use devices may create a risk of contamination and/or cause patient infection or cross-infection, including, but not limited to, the transmission of infectious disease(s) from one patient to another. Contamination of the device may lead to injury, illness, or death of the patient.
- Catheter ablation procedures present the potential for significant x-ray exposure, which can result in acute radiation injury as well as increased risk for somatic and genetic effects, to both patients and laboratory staff due to the x-ray beam intensity and duration of the fluoroscopic imaging. Catheter ablation should only be performed after adequate attention has been given to the potential radiation exposure associated with the procedure, and steps taken to minimize this exposure. Careful consideration must therefore be given for the use of the device in pregnant women.
- Implantable pacemakers and Automatic Implantable Cardioverter / Defibrillator (AICD) may be adversely affected by radiofrequency current. It is important to have temporary external sources of pacing available during ablation and temporarily reprogram the pacing system to minimum output or OFF mode to minimize risk of inappropriate pacing. Exercise extreme caution during ablation when in close proximity to atrial or ventricular permanent leads; program the AICD to the OFF mode during the ablation procedure; and, perform complete implantable device analysis on all patients after ablation.
- Do not immerse proximal handle or cable connector in fluids; electrical performance could be affected.
- Do not use excessive force to advance or withdraw the catheter when resistance is encountered.
- Do not scrub or twist the distal tip electrode during cleaning.
- The FireMagic™ Cardiac RF Ablation Catheter should be operated by appropriately trained personnel in a fully-equipped electrophysiology laboratory.
- Do not exceed 70 watts during ablation. If using a 100 watt generator, program the power to cut off at 70 watts.
- This device has not been shown to be safe at temperatures above 60 °C.
- Check connection reliability of accessory cable when weak signal or serious interference is encountered and replace a new one immediately when problem cannot be solved.
- When the catheter is connected to the supply mains-operated stimulator, inadvertently introduced leakage current may enter into the heart. Please use the internal power -operated simulator.
- Make sure that the output voltage used by the device is below rated accessory voltage of the catheter and its connection cable.
- The rated accessory voltage of catheter and Connection cable is 450V.
- When handling indwelling leads, the terminal pins or exposed metal are not to be touched nor be allowed to contact electrically conductive or wet surfaces.
- The device is disposable. Otherwise, MicroPort EP can not assure and control the occurrence of the potential risk due to the reuse of the device, like infection, loss of performance, loss of mechanical integrity.
- When used, catheter shall be disposed of as medical waste according local laws and regulations.

【PRECAUTIONS】

- The long-term risks of protracted fluoroscopy and creation of radiofrequency induced lesions have not been established. Careful consideration must therefore be given to the use of the device in prepubescent children. Furthermore, the risk/benefit in asymptomatic patients has not been studied.
- To avoid thromboembolism, usually intravenous heparin infusion should be given to maintain an activated clotting time>300 s during ablation or anticoagulation medicine is taken for about 3 months after the procedure as instructed by the doctor. No consensus is available yet as to the necessity for short-time anticoagulation after ablation.
- Careful catheter manipulation must be performed in order to avoid cardiac damage, perforation, or tamponade.
- Always pull the thumbknob back to straighten the catheter tip before insertion or withdrawal of the catheter.
- Apparent low power output, high impedance reading or failure of the equipment to function correctly at normal settings may indicate faulty application of the dispersive electrode(s) or failure of an electrical lead. Do not increase power before checking for obvious defects or misapplication.
- The sterile packaging and catheter should be inspected prior to use.
- It is recommended that less than thirty applications of energy be performed per catheter in 1 minute.
- The FireMagic™ Cardiac RF Ablation Catheter used in conjunction with a radiofrequency generator is capable of delivering significant electrical power. Patient or operator injury can result from improper handling of the catheter and dispersive

electrode, particularly when operating the device. During energy delivery, the patient should not be allowed to come in contact with grounded metal surfaces. If during ablation, tip temperature does not rise, discontinue ablation immediately and replace catheter.

- Regularly inspect and test re-usable cables and accessories. Carefully read this instruction for use before operating the catheter. Pay attention to all warnings and precautions, otherwise, accident may occur.
- Electromagnetic interference (EMI) produced by FireMagic used in conjunction with a radiofrequency generator during normal operation may adversely affect the performance of other equipment.
- Do not use near MRI equipment since movement or heating of the catheter may occur and the image on the display may become distorted.

【SIDE EFFECTS】

Complications that may occur during use of this device include pleural effusion; vascular arteriovenous fistula; local bleeding, vascular bleeding; hemothorax, pneumothorax, hemopneumothorax or chest pain; radiation injury; deep venous thrombosis, pulmonary; embolism, damage to the coronary vasculature and valves, valve damage; transient cerebral ischemic attack; pericardial effusion, cardiac tamponade, atrioventricular block, stroke, myocardial infarction, arterial perforation; vagal syncope at femoral puncture, transient traumatic or radiofrequency-related second and transient sinus bradycardia; death; Pericarditis; systemic infection.

【DEVICES USED IN COMBINATION】

- Transseptal sheath (8 French(F)) matching the catheter
- Radiofrequency (RF) Generator
- Cable for cardiac RF ablation catheter

【RADIOFREQUENCY ABLATION】

For radiofrequency ablation, the catheter must be connected to the appropriate input connector(s) on the radiofrequency generator. For proper RF generator interface, only use the connection cable by MicroPort EP.

NOTE: The catheter tip electrode must be switched from the electrogram recording equipment (via the generator controls or an external switch) to the radiofrequency generator power output for ablation (see applicable radiofrequency generator manual for instructions). For the temperature sensing catheter, verify that the generator displays a temperature near 37 °C after the catheter is inserted in the patient and before applying RF power.

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The FireMagic™ Cardiac RF Ablation Catheter is compatible with RF generators including HL-100F, IBI 1500T, Stockert EP Shuttle and Atakr II. Choose an appropriate connection cable of MicroPort in Table 2, and make the connection of RF generator and connection cable following the instructions in Table 3. Connection and operation instructions may also be found in the applicable connection cable manual. Refer to the applicable RF generator manual for detailed instructions as to generator operation.

Ask Shanghai MicroPort EP MedTech Co., Ltd. for more information about accessory cable types compatible with other RF generators.

NOTE: If radiofrequency current is interrupted for either a temperature or an impedance rise (a preset limit is exceeded), the catheter should be removed, and the tip cleaned of coagulum. When cleaning the tip electrode, be careful not to twist the tip electrode with respect to the catheter shaft; twisting may damage the tip electrode bond and loosen the tip electrode.

Table 2: Specification of connection cable

No	Cable Type	RF Generator	Connector (for RF Generator)
1	EPW B4	IBI1500T	
2	EPW ST	Stockert EP Shuttle, HL-100F	
3	EPW AK	Atakr II	

NOTE: 1. For details please refer to the applicable connection cable manual instructions.

2. The MAX output frequency for compatible RF generator should be at the range of 450KHZ-575KHZ.

Table 3: RF generators compatible with The FireMagic and its connection cable

HL-100F RF generator	
Direction for use	1. Connect The FireMagic™ Cardiac RF Ablation Catheter to HL-100F RF generator interface with appropriate connection cable; 2. Refer to the applicable HL-100F manual for detailed instructions as to generator operation for RF ablation.
Statement	This appendix is only applicable to the FireMagic™ Cardiac RF Ablation Catheter and the HL-100F RF generator used with it.
IBI 1500T generator	
Direction for use	1. Connect FireMagic™ Cardiac RF Ablation Catheter to IBI 1500T RF generator interface with appropriate connection cable of MicroPort;

	2. Refer to the applicable IBI 1500T manual for detailed instructions as to generator operation for RF ablation.
Statement	This appendix is only applicable to the FireMagic™ Cardiac RF Ablation Catheter and the IBI 1500T RF generator used with it.
Stockert EP Shuttle RF generator interface	
Direction for use	1. Connect the FireMagic™ Cardiac RF Ablation Catheter to Stockert EP Shuttle RF generator interface with appropriate connection cable. 2. Refer to the applicable Stockert EP Shuttle manual for detailed instructions as to generator operation for RF ablation.
Statement	This appendix is only applicable to the FireMagic™ Cardiac RF Ablation Catheter and the Stockert EP Shuttle RF generator used with it.
Atakr II RF generator	
Direction for use	1. Connect the FireMagic™ Cardiac RF Ablation Catheter to Atakr II RF generator interface with appropriate connection cable ; 2. Refer to the applicable Atakr II manual for detailed instructions as to generator operation for RF ablation.
Statement	This appendix is only applicable to the FireMagic™ Cardiac RF Ablation Catheter and the Atakr II RF generator used with it.

Note:

Please choose the neutral electrode which is recommended by the manufacturers of the compatible RF Generator, for the details, refer to the IFU of the generator, if not, please choose the neutral electrode obtained the CE approve.

【SUGGESTED INSTRUCTIONS FOR USE】

- (1) Choose appropriate catheter curve according to desired site within heart.
- (2) Check the catheter package and make sure there is no damage.
- (3) Remove catheter from package, carefully check the integrity of electrodes and catheter, and place in a sterile work area.
- (4) Create a vascular access in a large central vessel by using aseptic techniques and insert the catheter.
- (5) Connect the catheter to the recording equipment and/or the RF generator using the appropriate accessory cables.
- (6) Advance the catheter to the area of the endocardium which is to be examined. Use both fluoroscopy and electrograms to aid in proper positioning.
- (7) The catheter tip can be deflected to facilitate positioning by manipulating the thumbknob to vary tip curvature. Pushing the thumbknob forward causes the catheter tip to bend; when the knob is pulled back, the tip straightens. When it has been determined that the tip electrode is in stable contact with the intended ablation site, the catheter tip electrode connection must be switched from the recording equipment to the radiofrequency generator in preparation for delivery of radiofrequency current. Circuit impedance should be approximately 100 ohms upon initiation of radiofrequency current.
- (8) Radiofrequency current may be re-applied to the same or alternate sites using the same catheter. However, in the event of a generator cutoff (impedance or temperature), the catheter must be withdrawn and cleaned the coagulum on the tip electrode before radiofrequency current is re-applied. When the catheter is withdrawn, the tip electrode should be cleaned prior to re-insertion. A sterile gauze pad dampened with sterile saline may be used to gently wipe the tip section clean; do not scrub or twist the tip electrode, or damage to the tip electrode bond may occur and loosen the tip electrode.
- (9) If there are any questions regarding the use or performance of this product, please consult with the local distributor or the manufacturer.

【TRANSPORTATION REQUIREMENTS】

In transit, the product shall be protected from heavy load, direct sunlight and rain or as specified in the ordering contract. The temperature during transportation shall be kept between 0℃ and 45℃.

【STORAGE REQUIREMENTS】

The product shall be stored in a shady and cool, dry, clean, and well-ventilated warehouse which is in natural air circulation environment. The temperature during storage shall be kept between 0℃ and 45℃.

【SHELF LIFE】

The shelf life of product is three years while it meets the conditions for storage.

【STERILIZATION】

This product has been sterilized with ethylene oxide gas. Never re-sterilize and reuse it. Do not use the catheter if the package is open or damaged. Use the catheter prior to the expiration date shown on the package label.

【SYMBOLS EXPLANATION】

1.

**DO NOT REUSE**

2.

**CONSULT INSTRUCTIONS FOR USE**

3.

**PROTECT FROM HEAT SOURCE AND RADIATION SOURCE**

4.		KEEP DRY
5.		USE BY
6.		BATCH CODE
7.		CATALOGUE NUMBER
8.		STERILIZED USING ETHYLENE OXIDE
9.		DO NOT USE IF PACKAGE IS DAMAGED
10.		QUANTITY OF PRODUCT CONTAINED 1
11.		TEMPERATURE LIMITATION
12.		TYPE CF APPLIED PART
13.		DATE OF MANUFACTURE
14.		MANUFACTURER
15.		AUTHORISED REPRESENTATIVE IN THE EUROPEAN COMMUNITY
16.		IMPORTER

【AFTER-SALES SERVICE】

With “providing the medical sector with high quality and efficacy medical products” as its top operational objective, Shanghai MicroPort EP MedTech Co., Ltd. (hereinafter referred to as MicroPort EP Co.) guarantees that its products are free of defects in materials or manufacturing when the clients receive them. For other questions relating to the products, directly consult the company.

【SOLEMN STATEMENT】

Shanghai MicroPort EP MedTech Co., Ltd. expressly states herein that its **FireMagic™ Cardiac RF Ablation Catheter** is disposable and cannot be reused. MicroPort EP Co. will not recommend, indicate and imply in any manner the reusability of the system, and will not assume the responsibility for any accident or product damage resulting from reuse.

FireMagic™ Cardiac RF Ablation Catheter can be connected and used only with the compatible devices specified herein, and MicroPort EP Co. will not assume the responsibility for damage to product device, procedure failure and the like resulting from operating mistakes or any other human factors.

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Shanghai MicroPort EP MedTech Co., Ltd.

Address: Building 23&26&28, Lane 588, Tianxiong Rd., 201318 Shanghai, PEOPLE’S REPUBLIC OF CHINA

Postcode: 201318

Tel: +86 (21) 60969600

Fax: +86 (21) 20903925

E-mail: customerservice@everpace.com

Website: www.everpace.com

【AUTHORIZED REPRESENTATIVE IN THE EUROPEAN COMMUNITY】

MicroPort Medical B.V.

Address: Paasheuveweg 25, 1105BP Amsterdam, The Netherlands

E-mail: cs@microport-int.com

Tel: +31(0)205450100-8

Fax: +31(0)205450109

The copyright of this document is owned by Shanghai MicroPort EP MedTech Co., Ltd.
Copying is not allowed without its consent.

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**De copyright van dit document is eigendom van Shanghai MicroPort EP MedTech Co., Ltd.