Columbus™

使用说明书	中文	页码	1
体表参考电极			
Instructions for Use	English	Page	5
Columbus™ Reference Patch			
Kullanım Talimatları	Türkçe	Sayfa	8
Columbus™ Referans Bağlantısı			
Οδηγίες χρήσης	Ελληνικά	Σελίδα	11
Columbus™ επίθεμα αναφοράς			
Instrucciones de uso	Español	paginación	14
Parche de referencia Columbus™			
Instruções de Utilização	Português	Página	17
Emplastro de Referência Columbus™			
Инструкция за употреба	Български	Страница	20
Columbus™ референтна лепенка			
Инструкция по применению	Русский	Страница	23
эталонная накладка Columbus™			
Gebrauchsanweisung	Deutsch	Seite	26
Columbus™ Referenz-Patch			
Instructions d'utilisation	Français	Page	29
Patch de référence Columbus™			
Istruzioni per l'uso	Italiano	Pagina	32
Columbus™ Patch di Riferimento			
Instructies voor het gebruik	Nederlands	Pagina	35
Columbus™ referentiepatch			
Інструкція з застосування	Українська	Сторінка	38
Columbus™ Еталонна накладка			

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Shanghai MicroPort EP MedTech Co., Ltd.

【PRODUCT DESCRIPTION】

Columbus™ Reference Patch is an integral part of Columbus™ 3D EP Navigation System (hereinafter referred to as Columbus 3D) and Firemagic™ 3D Irrigated Ablation Catheter (hereinafter referred to as Firemagic 3D) for accurate locating and 3D graphics imaging.

The magnetic location sensors embedded in the adhesive patch provide the location information to the Columbus 3D by the FORLNK™ connection cable. The reference point which obtains by the patch can compensate and calibrate the accuracy about the location information by the 3D Irrigated Ablation Catheter.

【PRODUCT SPECIFICATION】

Table 1: Columbus™ Reference Patch Specification

No.	Model	Effective length
1	EPR 105 A	105 cm

【INDICATIONS】

The intended use of the Columbus™ reference patch is to provide a reference location relative to the 3D Irrigated Ablation Catheter when used in conjunction with the Columbus™ 3D EP Navigation System.

【CONTRAINDICATIONS】

None

【WARNING】

- Connection and use of this product should be performed by qualified physicians who have received professional training.
- Never immerse the near-end cable connector in liquid. Otherwise, the electric performance may be affected.
- When location is not displayed accurately or interference is increase, check interface reliability of connection cable. If abnormality cannot be excluded, the external reference patch should be replaced.
- When placing the catheter and connection cable, make sure that they are not in contact with the patient or other cables. When they are not in service for the time being, keep them in a place isolated from the patient.
- When used, catheter shall be disposed of as medical waste according local laws and regulations.
- The device is disposable. Otherwise, MicroPort EP can not assure and control the occurrence of the potential risk due to the reuse of the device, like infection, loss of performance, loss of mechanical integrity.

【PRECAUTIONS】

- Prior to use, carefully read the operating instructions for the product. Pay attention to all warnings and precautions in the instructions. Otherwise, accident may occur.
- Use with caution in patients with skin allergies.
- Do not paste the external reference patches directly to unhealed wounds.
- This product is only for a single-use. Do not re-sterilize and re-use.
- Before use, carefully observe the shelf life on package, do not use product with damaged packaging or expired product.

【SUGGESTED INSTRUCTIONS FOR USE】

- 1) Take the external reference patch out of package.
- 2) Tear off the backing on end electrode patch, paste the patch on heart region of patient's back with catheter body of external reference patch placed along the body long axis. Patch and catheter body can be steadily pasted to the back by medical tape, in order to avoid displacement by external force or patient's own movement.
- 3) Connect the product to Columbus 3D by EPS 150 C model connecting connection cable produced by the company.
- 4) Remove the product after operation and destroy according to local regulation.
- 5) Regarding any problem related to usage or product performance when using the product, please contact local distributor or manufacturer for consulting.

【TRANSPORTATION REQUIREMENTS】

In transit, the product shall be protected from heavy load, direct sunlight and rain or as specified in the ordering contract. The temperature during transportation shall be below 45℃.

【STORAGE REQUIREMENTS】

The product shall be stored in a shady and cool, dry, clean, and well-ventilated environment free of corrosive gases. The temperature during storage shall be kept in between 0℃ and 45℃.

【SHELF LIFE】

While meeting the conditions for storage, the shelf life is two years.

【STERILIZATION】

This product has been sterilized with ethylene oxide gas. Never re-sterilize and reuse it. Do not use the product if the package is open or damaged. Use the product prior to the expiration date shown on the package label.

【SYMBOLS EXPLANATION】

1.		DO NOT REUSE
2.		CONSULT INSTRUCTIONS FOR USE
3.		PROTECT FROM HEAT SOURCE AND RADIATION SOURCE
4.		KEEP DRY
5.		USE BY
6.		BATCH CODE
7.		CATALOGUE NUMBER
8.		STERILIZED USING ETHYLENE OXIDE
9.		DO NOT USE IF PACKAGE IS DAMAGED
10.		QUANTITY OF PRODUCT CONTAINED 1
11.		TEMPERATURE LIMITATION
12.		TYPE CF APPLIED PART
13.		DATE OF MANUFACTURE
14.		MANUFACTURER
15.		AUTHORISED REPRESENTATIVE IN THE EUROPEAN COMMUNITY
16.		IMPORTER

【AFTER-SALES SERVICE】

With “providing the medical sector with high quality and efficacy medical products” as its top operational objective, Shanghai MicroPort EP MedTech Co., Ltd. (hereinafter referred to as MicroPort EP Co.) guarantees that its products are free of defects in materials or manufacturing when the clients receive them. For other questions relating to the products, directly consult the company.

【SOLEMN STATEMENT】

Shanghai MicroPort EP MedTech Co., Ltd. (hereinafter referred to as MicroPort EP Co.) expressly states herein that its **Columbus™ Reference Patch** is disposable and cannot be reused. MicroPort EP Co. will not recommend, indicate and imply in any manner the reusability of the system, and will not assume the responsibility for any accident or product damage resulting from reuse.

Columbus™ Reference Patch can be connected and used only with the compatible devices specified herein, and MicroPort EP Co. will not assume the responsibility for damage to product device, procedure failure and the like resulting from operating mistakes or any other human factors.

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Shanghai MicroPort EP MedTech Co., Ltd.

Address: Building 23&26&28, Lane 588, Tianxiong Rd., 201318 Shanghai, PEOPLE'S REPUBLIC OF CHINA

Postcode: 201318

Tel: +86 (21) 60969600

Fax: +86 (21) 20903925

E-mail: customerservice@everpace.com

Website: www.everpace.com

【AUTHORIZED REPRESENTATIVE IN THE EUROPEAN COMMUNITY】

MicroPort Medical B.V.

Address: Paasheuvelweg 25, 1105BP Amsterdam, The Netherlands

E-mail: cs@microport-int.com

Tel: +31(0)205450100-8

Fax: +31(0)205450109

**The copyright of this document is owned by Shanghai MicroPort EP MedTech Co., Ltd.
Copying is not allowed without its consent.**

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**De copyright van dit document is eigendom van Shanghai MicroPort EP MedTech Co., Ltd.