FireMagic™ TrueForce™

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
Fig. 1-1: Normal type


Fig. 1-2: Micro-electrodes type

01-Tip Electrode; 02-Ring Electrodes; 03-Contact force sensor; 04-Magnetic field sensor; 05-Single lumen tube; 06-Deflectable section; 07- Magnetic field sensor; 08-Main body section; 9-Color code; 10-Thumbknob; 11-Handle; 12-PCB; 13-Connector; 14-Standard Luer; 15-Micro-electrodes

Fig. 1: The drawing of FireMagic™ TrueForce™ Ablation Catheter

FireMagic™ TrueForce™ Ablation Catheter is a multi-electrode luminal catheter with a deflectable tip designed to facilitate electrophysiological mapping of the heart and to transmit radiofrequency (RF) current to the catheter tip electrode for ablation purposes. The diameter of the catheter shaft is 8F. For ablation, the catheter is used in conjunction with an RF generator and a dispersive pad (indifferent electrode). The catheter has force-sensing technology that provides a real-time measurement of contact force between the catheter tip and the heart wall. The catheter has a high-torque shaft with a uni-directional deflectable tip section containing an array of electrodes which includes a 3.5 mm tip dome. All of the electrodes may be used for recording purposes. The tip electrode serves to deliver RF current from the RF generator to the desired ablation site. The tip electrode and ring electrodes are made from a noble metal. The catheter incorporates a thermocouple temperature sensor that is embedded in the 3.5 mm tip electrode. Tip deflection is controlled at the proximal end by a handpiece in which a piston slides; a thumbknob on the piston controls piston travel. When the thumbknob is pushed forward, the tip is deflected (curved). When the thumbknob is pulled back, the tip straightens. Currently, the available curves for the uni-directional version of the FireMagic™ TrueForce™ Ablation Catheter include B, D, F and J. At the proximal end of the catheter, a saline input port with a standard Luer fitting terminates from the open lumen. This saline port serves to permit the injection of normal saline to irrigate the tip electrode. During ablation, normal saline is passed through the internal lumen of the catheter and through the tip electrode, to irrigate and cool the ablation site as well as the electrode tip. An irrigation pump should be used to control the saline irrigation. The catheter interfaces with standard recording equipment and a compatible RF generator via accessory extension cables with the appropriate connectors. This catheter features a location sensor embedded in the tip section that transmits location information to Columbus™ 3D EP Navigation System.

In addition, there are four micro-electrodes in the tip electrode of micro-electrode type, which have the function of temperature measurement and ECG signal mapping to the tip electrode.

【PRODUCT SPECIFICATION】

FireMagic™ TrueForce™ Ablation Catheter has a color code near the handpiece. Different code colors show different shapes of the curve. For specifications and dimensions of the product available, see Table 1. Doctors can make choice in accordance with the location of disease and their own personal habits.

Table 1: Specification for FireMagic™ TrueForce™ Ablation Catheter

Model	Curve shape	Irrigati on type	Number of magnetic location sensor	Diamet er (mm)	Length (cm)	Micro-electrod es or not				
EPV8BTC	B/Red	6 holes	2	2.67	105	No				
EPV8DTC	D/Blue									
EPV8FTC	F/Orange									
EPV8JTC	J/Black									
EPV8BTM	B/Red	66 holes					2			
EPV8DTM	D/Blue									
EPV8FTM	F/Orange									
EPV8JTM	J/Black									
EPV8BTOC	B/Red	6 holes	1				No			
EPV8DTOC	D/Blue									
EPV8FTOC	F/Orange									
EPV8JTOC	J/Black									
EPV8BTOM	B/Red	66 holes	1							
EPV8DTOM	D/Blue									
EPV8FTOM	F/Orange									
EPV8JTOM	J/Black									
EPV8BZC	B/red	6 holes	2			Yes				
EPV8DZC	D/blue									
EPV8FZC	F/orange									
EPV8JZC	J/black									
EPV8BZM	B/red	66 holes	2							
EPV8DZM	D/blue									
EPV8FZM	F/orange									
EPV8JZM	J/black									
EPV8BZOC	B/red	6 holes	1				Yes			
EPV8DZOC	D/blue									
EPV8FZOC	F/orange									
EPV8JZOC	J/black									
EPV8BZOM	B/red	66 holes	1							
EPV8DZOM	D/blue									
EPV8FZOM	F/orange									
EPV8JZOM	J/black									

【INTENDED USE】

FireMagic™ TrueForce™ Ablation Catheter is indicated for use in cardiac electrophysiological mapping and, when used in conjunction with a radiofrequency generator, for cardiac ablation. FireMagic™ TrueForce™ Ablation Catheter provides a real-time measurement of contact force between the catheter tip and heart wall, as well as location information when used with Columbus™ 3D EP Navigation System.

【INTENDED CLINICAL BENEFITS】

When compared with 100% acute procedural success rate and 65% 12-month procedural success rate of non-CF-sensing catheter, CF-sensing catheter for the treatment of atrial fibrillation is safe and effective with 100% acute procedural success rate and 81.25% 12-month procedural success rate, which is clinically

meaningful improvement in acute outcome and long-term outcome, and a trend towards increased efficacy and a reduced complication rate.

【INDICATIONS】

Atrial fibrillation

【CONTRAINDICATIONS】

- Left atrium size $\geq 50\text{mm}$, or $\text{LVEF} \leq 40\%$;
- Left atrium thrombus;
- Patients once received atrioseptoplasty or had atrial myxoma;
- Those with cerebrovascular disease (including stroke and transient ischemic attack) in the past 3 months;
- The patient with cardiovascular events in the last 3 months (including acute myocardial infarction, coronary intervention or heart bypass surgery, prosthetic valve replacement or repair or atrial or ventricular incision);
- Patients with acute or severe systemic infection;
- Patients with severe combined liver, kidney diseases;
- Patients with obvious bleeding tendency and blood diseases;
- Patients with malignant tumors or end-stage diseases.

【TARGET GROUP】

For 18-75 years old adults, men or non-pregnant women.

【WARNING】

- With significant potential X-ray exposure, catheter ablation may subject patients and operators to acute radiation damage and increase the risk of impact on their bodies and genes. Catheter ablation can be applied only when sufficient attention is paid to the potential exposure involved in the procedure and measures are taken to minimize radiation exposure.
- Implantable pacemakers and implantable cardioverter/defibrillator (ICDs) may be adversely affected by RF current. It is important to: have temporary external sources of pacing and defibrillation available during ablation; temporarily reprogram the pacing system to minimum output or OFF mode to minimize the risk of inappropriate pacing. Exercise extreme caution during ablation when in close proximity to atrial or ventricular permanent leads; program the ICD to the OFF mode during the ablation procedure; and, perform complete implantable devices analysis on all patients after ablation.
- The temperature sensor in the catheter electrode is to measure the electrode temperature instead of the tissue temperature. Because of the cooling effect that normal saline infusion has on the electrode and the feature that RF ablation is passively heated, the electrode temperature displayed on RF generator is noticeably lower than that of the tissue in the ablation zone during the procedure. Temperature sensor is used to check if the infusion flow is appropriate, and at the same ablation temperature, larger infusion flow corresponds to higher RF discharge power and generates a deeper ablation zone. Therefore, special attention shall be paid to the infusion flow and actual power of ablation in procedure.
- Throughout the procedure, balance of patient's body fluid needs to be monitored to prevent overload of body fluid. Some patients, whose capability to handle overload of body fluid drops for some reasons, tend to suffer from acute pulmonary edema or heart failure during or after the procedure. This is particularly true of patients with congestive heart failure and elderly people. It is important to evaluate the level of risk that patients may develop overload of body fluid.
- Performing procedure again with this catheter on patients who has been subjected to RF ablation may face higher risk of perforation or pericardial effusion.
- When using Ablation Catheter, ferromagnetic devices and equipment will have potential offset interference with the product in picking up the signal for displaying ablation electrode and curve position. Therefore, throughout the procedure from point acquisition for 3D reconstruction to completion of ablation, approaching to or moving ferromagnetic devices and equipment such as MRI equipment and X-ray apparatus during procedure shall be avoided. In case an offset is present, usually returning the device to the original location will remove the offset.
- Never immerse the near-end handle or cable connector in liquid. Otherwise, the electric performance may be affected.
- This product shall be operated by trained personnel in electrophysiologically well-equipped catheter room.
- In insertion, never exert excessive force, and pull out the catheter in case of resistance.
- In cleaning, never wipe or twist the ablation end electrode heavily.
- To ensure patient's safety and the integrity of electrode catheter, never wipe the catheter with alcohol.
- Make sure that the output voltage used by the device is below rated accessory voltage of the catheter and connection cable.
- When placing the catheter and connection cable, make sure that they are not in contact with the patient or other cables. When they are not in service for the time being, keep them in a place isolated from the patient.
- When used, catheter shall be disposed of as medical waste according local laws and regulations.
- The catheter shall be connected to the non-implantable pulse generator before the pacing leads are connected to the catheter.
- When handling indwelling leads, the terminal pins or exposed metal are not to be touched nor be allowed to contact electrically conductive or wet surfaces.
- This device is packaged for single use only. Do not reuse, reprocess or resterilize. Reuse, reprocessing or resterilization may compromise the structural integrity of the device and/or lead to device failure that in turn may result in patient injury, illness or death. Also, reprocessing or resterilization of single use devices may create a risk of contamination and/or cause patient infection or cross-infection, including, but not limited to, the transmission of infectious disease(s) from one patient to another. Contamination of the device may lead to injury, illness, or death of the patient.
- If the patient with a history of heparin induced thrombocytopenia (HIT), bivalirudin is a preferred alternative to heparin for anticoagulation during procedures.
- When used with HF surgical equipment, parts of the equipment are provided with protection against burning of the patient.
- The contact force reading is for information only and is not intended to replace standard handling precautions.

- The catheter must be warmed up as specified prior to use. If the catheter has not reached a steady state condition, there is potential for a zero-offset drift to occur which could result in an inaccurate contact force reading.
- Always zero the contact force reading following insertion into the patient or when moving the catheter from one chamber of the heart to another. Ensure the catheter is not in contact with heart tissue prior to zeroing. Refer to the User Manual for your Columbus™ 3D EP Navigation System for instructions on how to zero the contact force reading.

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FireMagic™ TrueForce™ Ablation Catheter requires special precautions in electromagnetic compatibility (EMC), and must be installed and used in accordance with the EMC information provided in this Description. Portable and mobile RF communication devices will affect the operations of FireMagic™ TrueForce™ Ablation Catheter.

 **Warning:** Using unprescribed accessories, converters or cables other than the converters and cables which are approved as internal replacement parts by FireMagic™ TrueForce™ Ablation Catheter’s manufacturer, may lead to increased radiation or decreased anti-interference ability of the diagnostic catheter.

 **Warning:** FireMagic™ TrueForce™ Ablation Catheter should not be placed close to other equipment or be stacked. If the conditions above are not satisfied, you should observe and verify that it works normally with its working configuration.

Electromagnetic radiation						
FireMagic™ TrueForce™ Ablation Catheter is intended for use in the following specified electromagnetic environment. The customer or user should ensure that it works in such an environment.						
Radiation test	Compliance	Electromagnetic environment guidance				
Radiated emission CISPR 11: 2016	Group 1 Class A	FireMagic™ TrueForce™ Ablation Catheter its RF energy only for the internal functions in standby tests. Therefore, its RF emission is very low with a very small chance of interference from the nearby electronic equipment.				
Conducted emission CISPR 11: 2016	Not applicable					
Radiated emission CISPR 11: 2016	Group 1 Class A	FireMagic™ TrueForce™ Ablation Catheter is used in non-domestic equipment and facilities that are not directly connected to low-voltage public power grid in home dwellings.				
Voltage fluctuation / Scintillation emission IEC 61000-3-3:2017	Not applicable					
Harmonic radiation IEC 61000-3-2:2018	Not applicable					
Electromagnetic anti-interference						
FireMagic™ TrueForce™ Ablation Catheter is intended for use in the following specified electromagnetic environment. The customer or user should ensure that it works in such an environment.						
Anti-interference test	Test level	Working level				
Electrostatic discharge (ESD) IEC 60601-1-2:2014	±8kV contact discharge ± 15 kV air discharge	±8kV contact discharge ± 15 kV air discharge				
EFT IEC 60601-1-2:2014	±1kv 100kHz repetition frequency	±1kv 100kHz repetition frequency				
Surge IEC 60601-1-2:2014	Line-to-ground: ±2kv	Line-to-ground: ±2kv				
Voltage drop IEC 60601-1-2:2014	0% U _T ; 0.5 cycle: At 0°, 45°, 90°, 135°, 180°, 225°, 270°and 315° 0% U _T ; 1 cycle and 70% U _T ; 25/30 cycle: Single phase: at 0°	0% U _T ; 0.5 cycle: At 0°, 45°, 90°, 135°, 180°, 225°, 270°and 315° 0% U _T ; 1 cycle and 70% U _T ; 25/30 cycle: Single phase: at 0°				
Voltage Interruptions	0% U _T ; 250/300 cycle	0% U _T ; 250/300 cycle				
Power frequency(50 Hz) Magnetic field IEC 60601-1-2:2014	30A/m	30A/m				
Note: UT is the AC power supply voltage before applying the corresponding test level. Test FireMagic™ TrueForce™ Ablation Catheter in 100 and 230 VAC.						
Electromagnetic anti-interference						
FireMagic™ TrueForce™ Ablation Catheter is intended for use in the following specified electromagnetic environment. The customer or user of the deflectable curve mapping catheter should ensure that it works in such an environment.						
Anti-interference Test	Test Level	Compliance Level	Electromagnetic environment Guidance			
Conductive RF IEC 60601-1-2:2014	0.15MHz-80MHz	3 V	The portable and mobile RF communication devices shall not be used within a certain distance to any part of FireMagic™ TrueForce™ Ablation Catheter, including the cable and the recommended distance may be calculated by the equation applicable to the transmitter frequency.			
	In ISM bands between 0.15MHz and 80MHz 80% AM at 1kHz	6V				
Radiation RF IEC 60601-1-2:2014	80-2.7GHz	3 V/m				
Test specifications for ENCLOSURE PORT IMMUNITY to RF wireless communications equipment						
Test frequency (MHz)	Band ^{a)} (MHz)	Service ^{a)}	Modulation ^{b)}	Maximum power (W)	Distance (m)	IMMUNITY TEST LEVEL
385	280-390	TETRA 400	Pulse modulation ^{b)}	1.8	0.3	27

			18Hz			
450	430-470	GMRS 460, FRS 460	FM ^{c)} ± kHz deviation 1kHz sine	2	0.3	28
710	704-787	LTE Band 13, 17	Pulse Modulation ^{b)} 217 Hz	0.2	0.3	9
745						
780						
810	800-960	GSM 800/900, TETRA 800, iDEN 820, CDMA 850, LTE Band 5	Pulse Modulation ^{b)} 18 Hz	2	0.3	28
870						
930						
1720	1700- 1990	GSM 1800; CDMA 1900; GSM 1900; DECT; LTE Band 1, 3, 4, 25; UMTS	Pulse Modulation ^{b)} 217 Hz	2	0.3	28
1845						
1970						
2450	2400- 2570	Bluetooth, WLAN, 802, 11 b/g/n, RFID 2450, LTE Band 7	Pulse Modulation ^{b)} 217 Hz	2	0.3	28
5240	5100- 5800	WLAN 802, 11 a/n	Pulse Modulation ^{b)} 217 Hz	0.2	0.3	9
5500						
5785						
Note: If necessary to achieve the IMMUNITY TEST LEVEL, the distance between the transmitting antenna and FireMagic™ TrueForce™ Ablation Catheter may be reduced to 1 m. The 1 m test distance is permitted by IEC 61000-4-3.						

【PRECAUTIONS】

- In case of obvious low power output, readings of high impedance, abnormal signal disturbing and the like, immediately check if the dispersive pad is stuck and damaged, the infusion tube unblocked, the external reference patch securely bonded, and the connecting cable reliably connected; and if the abnormality cannot be removed, do not increase the ablation power.
- To avoid thromboembolism, usually intravenous heparin infusion should be given to maintain an activated clotting time>300 s during ablation or anticoagulation medicine is taken for about 3 months after the procedure as instructed by the doctor, giving NOAC priority. No consensus is available yet as to the necessity for short-time anticoagulation after ablation.
- Flush the catheter and infusion tube. Before pushing the catheter into the body, perform complete degassing to prevent air embolism, and during the service, always maintain certain low infusion flow rate to prevent the blood from coagulation in the infusion tube.
- Make sure to carefully follow the power-adjusting steps specified in the operating instructions. Excessively increasing the ablation power may cause steam burst and perforation.
- While operating the catheter, always be careful to avoid damage to the heart which causes perforation or tamponade.
- Before pushing in or pulling out the catheter, always fully pull back the thumbknob so as to straighten the catheter.
- The catheter used in conjunction with RF ablation device can transfer considerable electric energy. Improperly operating the catheter and dispersing electrode (especially when operating the device) may damage the patient or operator. During the energy transfer, patient must not touch any earthed metal surface. If the tip temperature does not rise during ablation, immediately stop ablation and replace the catheter.
- Regularly check and test the reusable Connection cable and accessories.
- For connection, use equipment and cables specified in the operating instructions and compatible with this device.
- Prior to use, carefully read the operating instructions for the product. Pay attention to all warnings and precautions in the instructions. Otherwise, accident may occur.
- Prior to use, the sterile package and catheter shall be inspected. Remove the catheter from its package carefully and place it in a sterile work area.
- Do not use near MRI equipment since movement or heating of the catheter may occur and the image on the display may become distorted.
- Conductive parts of electrodes and associated connectors for applied parts including the neutral electrode not contact other conductive parts and earth.
- Do not use the catheter if the package is open or damaged. Use the catheter prior to the expiration date shown on the package label.

【ADVERSE REACTIONS】

Inconvenient use of the product, pain, discomfort, transient ischemic attack, inflammatory reaction, toxic reaction, infection, radiation injury, phrenic nerve palsy, cardiac perforation/cardiac tamponade, pericardial effusion, pulmonary vein stenosis, hemothorax, pneumothorax, thrombosis/aeroembolism, stroke/cerebrovascular event, pulmonary edema, pericarditis/endocarditis, bleeding/vascular complications(arteriovenous fistula, hematoma, pseudoaneurysm), stiff left atrial syndrome, myocardial infarction, atrial esophageal fistula, phrenic nerve injury, heart block, valve injury, death.

【DEVICES USED IN COMBINATION】

- Transseptal Sheath matching the catheter
- Radiofrequency (RF) Generator
- Cable for TrueForce™ Ablation Catheter

- Columbus™ 3D EP Navigation System
 - Irrigation Pump
 - Tubing Set
 - External reference patch
 - Connection Box
- 【CONNECTION FOR PROCEDURE】**


Fig. 2-1: The sketch for connection (for normal type)

4. Catheter 5. Catheter handle 6. Tubing Set 7. Standard Luer
8. Cable for TrueForce™ Ablation Catheter 9. Connection cable
10. Connection Cable 11. Columbus™ external reference patch
12. Connector for the external reference patch


Fig. 2-2: The sketch for connection (for micro-electrode type)

1. Catheter 1-2. Catheter handle 1-3. Standard Luer
2. Irrigation pump 2-1. Tubing Set
3. Connection box 3-1. Connector of the catheter 3-2. Connection for micro-electrodes
4. Connection cable for ECG
5. PIU 5-1. External reference patch connector 5-2. Connector for catheter
5-3. ECG connector 5-4. Connector for RF generator
6. Cable
7. External reference patch 7-1. Columbus™ external reference patch 7-2. Connector for the external reference patch
8. RF generator 8-1. RF interface 8-2. Multi-thermocouple interface
9. Cable for connecting PIU and RF generator

As shown in Fig. 2, connect FireMagic™ TrueForce™ Ablation Catheter with relevant procedure equipment. Equipment compatible or used in combination include:

Compatible Equipment	Brand and Type	Remarks
RF Generator	OptimAblate™ EPE-RCAF-1A	It should be connected through Columbus™ 3D EP Navigation System
Irrigation Pump	OptimAblate™ EPE-IGP-1A	/
3D EP Navigation System	Columbus™ EPE-SYS-1A/EPE-SYS-2A	/
Reference Patch	Columbus™ EPR 105A	
External Reference Patch	Columbus™ EPTR150A, EPTR150B, EPTR150C, EPTR150D, EPR105B, EPRT75A	

Note: 1) The MAX output frequency for compatible RF generator should be at the range of 450KHZ-575KHZ

2) Please choose the neutral electrode which is recommended by the manufacturers of the compatible RF Generator, for the details, refer to the IFU of the generator, if not, please choose the neutral electrode obtained the CE approve.

• Connection Cable

The cable is chosen and used for connection of FireMagic™ TrueForce™ Ablation Catheter with 3D electrophysiological mapping system. Via connection cable, FireMagic™ TrueForce™ Ablation Catheter can interface with external equipment. The connection cable is the accessory for FireMagic™ TrueForce™ Ablation Catheter.

Name	Types	Remarks
Cable for TrueForce™ Ablation Catheter	EPV 150A	/
	EPV 180A	
	EPV 150B	
Connection box	EPBE-2A	For Micro-electrodes type

Note: The cable is the product of MicroPort EP.

The rated accessory voltage of catheter and connection cable is 250Vp

• Connection with Columbus™ 3D EP Navigation System

Select specified type of cable. Based on the corresponding “key” position, plug the plug on the side marked “catheter” of the connection cable in the jack in the catheter, and a “click” indicates the completion of the plugging; and plug the plug on the other side of the cable in the jack marked “Map interface” in the PIU of Columbus™ 3D EP Navigation System based also on the

corresponding “key” position and a “click” indicates the completion of the plugging.

Remove the cable from the connecting devices on the two sides: just pull it out with slight force by holding the arrow on the cable connector.

- **Connection with irrigation pump**

Butt and tighten the standard Luer on the compatible irrigation pump strap to the saline tube standard Luer on the handle of FireMagic™ TrueForce™ Ablation Catheter.

- **Connection with other devices and equipment**

1) For the connection and use of external reference patch, refer to IFU for Columbus™ External Reference Patch;

2) For the connection with and use of Columbus™ 3D EP Navigation System, refer to IFU for Columbus™ 3D EP Navigation System.

【SUGGESTED INSTRUCTIONS FOR USE】

1. To choose a catheter, note that the shape of the curve suits the heart anatomy to be reached;
2. First, carefully check if the package of the catheter is intact;
3. Remove the catheter from the package, carefully check the integrity of the electrodes and catheter, and place them on the sterile bench;
4. Establish channel in the peripheral great vessels by sterile technology;
5. In order to prevent damage to the catheter tip, use the insertion tube supplied with the catheter to advance or retract the catheter through the hemostasis valve of the sheath. After insertion, slide the insertion tube back toward the handle;
6. To verify compatibility between the sheath and catheter, advance the catheter through sheath prior to insertion;
7. Screw the catheter standard Luer to the strap connector of the infusion pump, vent by opening to the max infusion flow rate ($> 30\text{ml/min}$) and check that the infusion channel is unblocked: all the saline irrigation holes in the ablation electrode bleed water normally;
8. Switch the Irrigation pump flow rate to and maintain it at the low flow rate without ablation ($2\text{-}5\text{ ml/min}$), and insert the catheter into the established vessel channel;
9. Connect the catheter to the RF generator or 3D equipment with compatible connecting cable, and check that the connection with the RF generator, multi-channel electrophysiological recorder, stimulator and the like is properly performed;
10. Push the catheter into the endocardium zone to be examined and use X-ray and ECG to help with correct locating;
11. In order to achieve optimal force reading accuracy and stability, allow the catheter to warm up for 15 minutes after connection to the Columbus™ 3D EP Navigation System, prior to use of the force feedback feature. And allow the catheter to warm up for 5 minutes after inserted into the patient;
12. Zero the contact force reading following insertion into the patient. The tip electrode and first two ring electrodes must be outside of the sheath so that the force sensor is inside the body. Ensure the catheter tip is not in contact with tissue by evaluating the location on fluoroscopy and the Columbus™ System, the EGM amplitude, and catheter movement. Variations in the force reading at the same rate as the cardiac or respiration cycle may indicate contact with cardiac structures. Once these markers indicate the tip is not in contact, the reading can be zeroed. Refer to the User Manual for your Columbus™ 3D EP Navigation System for instructions on how to zero the contact force reading;
13. Zero the contact force reading when moving the catheter from one chamber of the heart to another or upon reinsertion;
14. Change the bending of the tip with the thumbknob to help locating by the catheter tip. When the thumbknob is pushed forward, the tip is deflected. When the thumbknob is pulled back, the tip straightens;
15. Operate the catheter so that the tip electrode arrives at the target location in the heart cavity, perform point acquisition for 3D reconstruction (when using the FireMagic™ and connecting the Columbus™ system) or pick up ECG;
When the position of the tip electrode or catheter bending cannot be displayed or the data is inaccurate, check if the connecting interfaces are correctly stabilized; or replace the cable or catheter. If the replacement results in returning to normal (the normal condition), it indicates the damage to the device;
16. After selecting the target location for ablation, set the irrigation pump flow rate and ablation parameters, switch to high irrigation flow rate and start ablation; Recommendation for selection of infusion flow: prior to ablation, infuse at high flow for 5s, and maintain the high flow within 5s at the completion of ablation; for an ablation power of 15-30W, use a high flow at 17ml/min , and for an ablation power of 30-50W, use a high flow at 30ml/min .
Before applying RF energy, verify that the infusion flow is increased by lowering the ablation electrode temperature;
And the duration for each ablation shall be within 120s. During ablation, it is allowed to move the ablation electrode to the next location;
When starting the RF current, the circuit impedance shall be about 100ohms; During discharge, carefully check the output power. An excessively low output power may suggest that the Saline irrigation holes in the ablation electrode are blocked. At that time, discharge must stop and the catheter must be pulled out. Clean the coagulum in the ablation electrode, verify that the Saline irrigation holes bleed water normally, and re-insert after complete venting is done: clean the tip gently with sterile gauze pad soaked with sterile saline; and never rub or twist the ablation electrode heavily. Otherwise, the bonding of the ablation electrode may be damaged, which leads to the electrode loosening. Alternatively, flush the infusion tube by injecting normal saline at high pressure with a syringe.
17. At the completion of the procedure, remove the catheter and dispose of it as specified.
18. In case of any questions on usage or product performance, consult local distributor or the manufacturer.

【TRANSPORTATION REQUIREMENTS】

In transit, the product shall be protected from heavy load, direct sunlight and rain or as specified in the ordering contract.

【STORAGE REQUIREMENTS】

The product shall be stored in a shady and cool, dry, clean, and well-ventilated warehouse which is in natural air circulation environment. The temperature during storage shall be kept between 0℃ and 45℃.

【SHELF LIFE】

The shelf life of product is three years while it meets the conditions for storage.

【STERILIZATION】

This product has been sterilized with ethylene oxide gas. Never re-sterilize and reuse it. Do not use the catheter if the package is open or damaged. Use the catheter prior to the expiration date shown on the package label.

【SYMBOLS EXPLANATION】

1.		DO NOT REUSE
2.		CONSULT INSTRUCTIONS FOR USE
3.		PROTECT FROM HEAT SOURCE AND RADIATION SOURCE
4.		KEEP DRY
5.		USE BY
6.		BATCH CODE
7.		CATALOGUE NUMBER
8.		STERILIZED USING ETHYLENE OXIDE
9.		DO NOT USE IF PACKAGE IS DAMAGED
10.		QUANTITY OF PRODUCT CONTAINED 1
11.		TEMPERATURE LIMITATION
12.		TYPE OF APPLIED PART
13.		DATE OF MANUFACTURE
14.		MANUFACTURER
15.		AUTHORISED REPRESENTATIVE IN THE EUROPEAN COMMUNITY
16.		UK RESPONSIBLE PERSON
17.		UNIQUE DEVICE IDENTIFIER
18.		IMPORTER

【SUPPLEMENTARY INFORMATION】

- Eudamed is available at: <https://ec.europa.eu/tools/eudamed>.
- The SSCP of FireMagic™ TrueForce™ Ablation Catheter is available at: <https://www.everpace.com>
- Electronic IFU is available at: <https://www.everpace.com>
- Basic UDI-DI: 69524509C6A04FL

【AFTER-SALES SERVICE】

With “providing the medical sector with high quality and efficacy medical products” as its top operational objective, Shanghai MicroPort EP MedTech Co., Ltd. (hereinafter referred to as MicroPort EP Co.) guarantees that its products are free of defects in materials or manufacturing when the clients receive them. For other questions relating to the products, directly consult the company.

【SOLEMN STATEMENT】

Shanghai MicroPort EP MedTech Co., Ltd. expressly states herein that its FireMagic™ TrueForce™ Ablation Catheter is disposable and cannot be reused. MicroPort EP Co. will not recommend, indicate and imply in any manner the reusability of the system, and will not assume the responsibility for any accident or product damage resulting from reuse.

FireMagic™ TrueForce™ Ablation Catheter can be connected and used only with the compatible devices specified herein, and MicroPort EP Co. will not assume the responsibility for damage to product device, procedure failure and the like resulting from operating mistakes or any other human factors.

Any serious incident that occurred in relation to the device should be reported to the manufacturer and the competent authority of the Member State.

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Shanghai MicroPort EP MedTech Co., Ltd.

Address: Building 23&26&28, Lane 588, Tianxiong Rd., 201318 Shanghai, PEOPLE'S REPUBLIC OF CHINA

Postcode: 201318

Tel: +86 (21) 60969600

Fax: +86 (21) 20903925

E-mail: customerservice@everpace.com

Website: www.everpace.com

【AUTHORIZED REPRESENTATIVE IN THE EUROPEAN COMMUNITY】

MicroPort Medical B.V.

Address: Paasheuvelweg 25, 1105BP Amsterdam, The Netherlands

E-mail: cs@microport-int.com

Tel: +31(0)205450100-8

Fax: +31(0)205450109

【UK RESPONSIBLE PERSON】

SUNGO Certification Company Limited

Address: 3rd floor, 70 Gracechurch Street, London, EC3V 0HR, England, United Kingdom

Contact person: Jackie Zhu (Ms.)

E-mail: ukrp@sungoglobal.com/Jackie.zhu20@outlook.com


Tel: +44-(0) 7868 234015

**The copyright of this document is owned by Shanghai MicroPort EP MedTech Co.,
Ltd.
Copying is not allowed without its consent.**

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
The diagram illustrates the setup of the EP navigation system for a patient. A patient is shown lying on an operating table (11). The system components and their connections are as follows:

- 1. Irrigatiepomp** (Irrigation pump) is connected to the catheter (6) via a tube (4).
- 2. 3D EP-navigatiesysteem** (3D EP navigation system) is connected to the catheter (6) via a tube (5) and to the RF-generator (3) via a cable (8).
- 3. RF-generator** (RF generator) is connected to the catheter (6) via a cable (9) and to the operating table (11) via a cable (10).
- The catheter (6) is inserted into the patient's heart (7) and connected to the RF-generator (3) via a cable (8).
- The operating table (11) is connected to the RF-generator (3) via a cable (10) and to the 3D EP navigation system (2) via a cable (12).

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De copyright van dit document is eigendom van Shanghai MicroPort EP MedTech Co., Ltd.