EasyFinder™

使用说明书 EasyFinder™ 一次性使用固定弯标测导管	中文	页码	1
Instructions for Use EasyFinder™ Electrophysiology Diagnostic Catheter	English	Page	11
Kullanım Talimatları EasyFinder™ Elektrofizyoloji Tanısal Kateteri	Türkçe	Sayfa	16
Οδηγίες χρήσης EasyFinder™ Διαγνωστικός καθετήρας ηλεκτροφυσιολογίας	Ελληνικά	Σελίδα	21
Instrucciones de uso Catéter de diagnóstico electrofisiológico EasyFinder™	Español	Página	26
Instruções de Utilização Cateter Diagnóstico de Eletrofisiologia EasyFinder™	Português	Página	31
Инструкция за употреба EasyFinder™ електрофизиологичен (ЕФ) диагностичен катетър	Български	Страница	36
Инструкция по применению Катетердиагностический электрофизиологический EasyFinder™	Русский	Страница	41
Bedienungshinweise EasyFinder™ Elektrophysiologie Diagnostikkatheter	Deutsch	Seite	46
Instructions d'utilisation Cathéter de diagnostic en électrophysiologie EasyFinder™	Français	Page	51
Istruzioni per l'uso EasyFinder™ catetere diagnostico per elettrofisiologia	Italiano	Pagina	56
Instructies voor het gebruik EasyFinder™ elektrofysiologische diagnostische katheter	Nederlands	Pagina	61
Інструкція з застосування EasyFinder™ Діагностичний катетер	Українська	Сторінка	66

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Shanghai MicroPort EP MedTech Co., Ltd.

【CATHETER DESCRIPTION】

EasyFinder™ Electrophysiology Diagnostic catheter (EasyFinder for short) manufactured by Shanghai MicroPort EP MedTech Co., Ltd. (MicroPort EP for short) is divided into two categories, one is the fixed curve diagnostic catheter and the other is the steerable curve diagnostic catheter.

EasyFinder has a high-torque shaft with an array of platinum electrodes at the distal tip. The platinum electrodes are used in cardiac electrophysiological mapping, recording and stimulation.

The tip section of the fixed curve diagnostic catheter has a preformed curve.


Fig. 1: EasyFinder™ Fixed Curve Diagnostic Catheter

The tip section of the steerable curve diagnostic catheter can be deflectable. The tip deflection is controlled at the proximal end by a thumbknob in which a piston sliders. When the piston is pushed forward with the thumbknob, the tip is curved. When the piston is pulled back, the tip straightens.


Fig. 2: EasyFinder™ Steerable Curve Diagnostic Catheter

At the proximal end of the catheter, a connector is provided for connecting specified cable. The connection cable as an important accessory of diagnostic catheter, is a electrical transmission line connecting catheter and standard recording equipment.

【PRODUCT SPECIFICATION】

The EP diagnostic catheter has many specifications for doctor to choose in accordance with the size of cardiac chambers, the approaching path and the site to be mapped.

Different specifications of catheters share the same principle, basic structure, materials, and intended use. Please refer to the following tables to get the detailed information for the specifications.

Table 1: Specifications of EasyFinder™ Fixed Curve Diagnostic Catheter

Specification	Curve/color	Electrode Num. & Electrodes spacing (MM)	Tip electrode	Length of shaft (CM)	Diameter of shaft (MM)	Recommend using site
EPD6AD252	A/Yellow	4 & 2-5-2	None	115	2.0	Right atrial (RA)
EPD6AD005		4 & 5-5-5				
EPD6AD010		4 & 10-10-10				
EPD6DD252	D/Blue	4 & 2-5-2				Right Ventricle (RV)
EPD6DD005		4 & 5-5-5				
EPD6DD010		4 & 10-10-10				
EPD6FD252	F/Black	4 & 2-5-2		60	1.65	His Bundle (HIS)
EPD6FD005		4 & 5-5-5				
EPD6FD010		4 & 10-10-10				
EPD5PE282	P/Grey	10 & 2-8 in sequence				Coronary sinus

English

EPD6PE282						(CS)
EPD6DE282	D/Blue			92		
EPD6AB252	A/Yellow	4 & 2-5-2	One	115	2	RA
EPD6AB005		4 & 5-5-5				
EPD6AB010		4 & 10-10-10				
EPD6DB252	4 & 2-5-2	RV				
EPD6DB005	4 & 5-5-5					
EPD6DB010	4 & 10-10-10					
EPD6FB252	4 & 2-5-2	HIS				
EPD6FB005	4 & 5-5-5					
EPD6FB010	4 & 10-10-10					
EPD5PC282	P/Grey	10 & 2-8 in sequence		60	1.65	CS
EPD6PC282	P/Grey	10 & 2-8 in sequence		60	2	
EPD6DC282	D/Blue	10 & 2-8 in sequence		92		
EPD5CG010	C/Green	2&10	None	115	1.65	RA
EPD5JG010	J/Orange	2&10			2	
EPD6CG010	C/Green	2&10				
EPD6JG010	J/Orange	2&10				

Table 2: Specifications of EasyFinder™ Steerable Curve Diagnostic Catheter

Specification	Curve/Color	Electrode Num. & Electrodes spacing (MM)	Tip electrode	Length (CM)	Diameter of shaft (MM)	Recommend using site
EPJ6DH252	D/Blue	4 & 2-5-2	None	115	2.0	RA/RV/HIS
EPJ6DH005		4& 5-5-5				
EPJ6DH010		4&10-10-10				
EPJ6DH282		10 & 2-8 in sequence				

Table 3: Corresponding table about the color code, shape of the curve

Shape of the curve	D	A	F	P	J	C
Codes color	Blue	Yellow	Black	Grey	Orange	Green
Schematic for the curve						

【INTENDED USE】

EasyFinder™ Electrophysiology Diagnostic Catheter is indicated for use in electrogram recording and cardiac stimulation during cardiac diagnostic electrophysiology studies for arrhythmia patients.

【INDICATIONS】

Cardiac diagnostic electrophysiology studies for arrhythmia patients

【CONTRAINDICATIONS】

- 1) Patients who have intracardiac mural thrombus or have had a ventriculotomy or atriotomy within the preceding four weeks.
- 2) Patients with prosthetic valves.
- 3) Active systemic infection.
- 4) Patients with obvious hemorrhagic tendency and have hematopathy.
- 5) Patients with malignant tumor and have terminal disease.
- 6) Patients with advanced heart failure (NYHA Functional Class III-IV).
- 7) Unstable angina and acute myocardial infarction within three months.
- 8) Stroke and transient ischemic attack within the last two weeks.

【TARGET GROUP】

For 18-80 years old adults, men or non-pregnant women

【WARNING】

- Do not attempt to operate the EasyFinder prior to completely reading and understanding these directions for use.
- This device is packaged for single use only. Do not reuse, reprocess or resterilize. Reuse, reprocessing or resterilization may compromise the structural integrity of the device and/or lead to device failure that in turn may result in patient injury, illness or death. Also, reprocessing or resterilization of single use devices may create a risk of contamination and/or cause patient infection or cross-infection, including, but not limited to, the transmission of infectious disease(s) from one patient to another. Contamination of the device may lead to injury, illness, or death of the patient.

- Do not use if the package is open or damaged; Use the device prior to the "Use By" data on the package label.
- Cardiac catheterization procedures present the potential for significant x-ray exposure, which can result in acute radiation injury as well as increased risk for somatic and genetic effects, to both patients and laboratory staff due to the x-ray beam intensity and duration of the fluoroscopic imaging. Cardiac catheterization should only be performed after adequate attention has been given to the potential radiation exposure associated with the procedure, and steps taken to minimize this exposure. Careful consideration must be given for the use of this catheter in pregnant women.
- Do not immerse the proximal handle or cable connector in fluids; electrical performance could be affected.
- Do not autoclave the catheter
- Always pull the thumb knob of the catheter back before insertion or withdrawal to assure that the catheter tip assumes its original shape.
- Careful catheter manipulation must be performed in order to avoid cardiac damage, perforation, or tamponade. Catheter advancement and placement should be done under fluoroscopic guidance. Do not use excessive force to advance or withdraw the catheter when resistance is encountered.
- To ensure patient safety and electrode catheter integrity, do not wipe the catheter with alcohol, or try to bend tube shaping.
- The conductive part of the catheter and cable should not be contact with other conductors, including earth.
- Cardiac catheterization procedures should be performed by appropriately trained personnel in a fully equipped electrophysiology laboratory.
- When used, catheter shall be disposed of as medical waste according to local laws and regulations.
- When the catheter is connected to the supply mains-operated stimulator, inadvertently introduced leakage current may enter into the heart. Please use the internal power -operated stimulator.
- Before handling the external pulse generator, the patient cable or indwelling leads, steps shall be taken to equalize the electrostatic potential between the user and the patient, for example by touching the patient at a site remote from the pacing lead.
- The catheter shall be connected to the non-implantable pulse generator before the pacing leads are connected to the catheter.
- When handling indwelling leads, the terminal pins or exposed metal are not to be touched nor be allowed to contact electrically conductive or wet surfaces.
- The device is disposable. Otherwise, MicroPort EP can not assure and control the occurrence of the potential risk due to the reuse of the device, like infection, loss of performance, loss of mechanical integrity.
- Do not use near MRI equipment since movement or heating of the catheter may occur and the image on the display may become distorted.

【ADVERSE REACTIONS】

Complications that may occur during use of this device include valve injury; vascular bleeding, local hematomas; thrombosis, thrombembolia; stroke; pulmonary embolism; pericarditis; vascular arteriovenous fistula, pseudoaneurysm; pneumothorax and hemothorax; myocardial infarction; cardiac tamponade; cardiac perforation; death; radiation injury; systemic infection.

【DEVICES USED IN COMBINATION】

- Transseptal sheath (6 French(F), 7 French(F)) matching the catheter
- Cable for mapping catheter



【CONNECTION CABLE】

The connection cable (Cable for mapping catheter) is an important accessory of the catheter.

It designed to connect the EP diagnostic catheter and the appropriate stimulation and/or standard recording equipment.

There are two cable types according to the pin number (see table 4). The pins (or sub-line seat) are marked with numbers D, 2, 3, and etc, which is corresponding to the electrode order on the EP diagnostic catheter.

Table 4: Specification of connection cable for EP diagnostic catheter

Specification	Number of Electrodes	Cable Parameter	
		Number of pins	Connector (for standard recording equipment)
COG010	10	10	
COG004	2 or 4	4	

Note: See the INSTRUCTIONS FOR USE for more detailed information about the connection cable which produced by MicroPort EP.

【SUGGESTED INSTRUCTIONS FOR USE】

- 1) Select the suitable catheter curvature according to the mapping parts of the heart;
- 2) Check carefully if there is any damage to the catheter package;
- 3) Take out the catheter from the package carefully without violent pulling out, for the fixed curve diagnostic catheter, take off the front curved-retaining tube, carefully check the integrity of the electrode and the catheter, and place it in a sterile work area;
- 4) For the steerable curve diagnostic catheter, please confirm that the thumbknob is pulled back completely before insertion. Advance the catheter to the area of the endocardium under investigation. Use both fluoroscopy and electrograms to aid in

- proper positioning. Adjust the radius of curvature as necessary by manipulating the thumbknob. Pushing the thumbknob forward causes the catheter tip to curve; when the knob is pulled back, the tip straightens;
- 5) Using puncture to build an access in the central large vessel, insert the catheter;
 - 6) Advance the catheter to the endocardial area to be examined. Please do it Under the X-ray and ECG in order to position the catheter correctly;
 - 7) Use the appropriate cable to connect the catheter with the standard recording device;
 - 8) Remove and dispose the catheter in an appropriate manner when the operation is finished. Do not re-sterilize and re-use. For the steerable curve diagnostic catheter, prior to removal of the catheter, confirm that the thumbknob has been pulled back completely. Remove the catheter and dispose of it in an appropriate manner. Do not resterilize and reuse;
 - 9) If there are any questions regarding the use of performance of this product, please consult with the local distributor of the manufacturer.

【TRANSPORTATION REQUIREMENTS】

In transit, the product shall be protected from heavy load, direct sunlight and rain or as specified in the ordering contract. The temperature during transportation shall be kept between 0°C and 45°C.

【STORAGE REQUIREMENTS】

The product shall be stored in a shady and cool, dry, clean, and well-ventilated warehouse which is in natural air circulation environment. The temperature during storage shall be kept between 0°C and 45°C.

【SHELF LIFE】

The shelf life of EasyFinder™ Fixed Curve Diagnostic Catheter is three years while it meets the conditions for storage.

The shelf life of EasyFinder™ Steerable Curve Diagnostic Catheter is two years while it meets the conditions for storage.

【STERILIZATION】

This product has been sterilized with ethylene oxide gas. Never re-sterilize and reuse it. Do not use the catheter if the package is open or damaged. Use the catheter prior to the expiration date shown on the package label.

【SYMBOLS EXPLANATION】

1.		DO NOT REUSE
2.		CONSULT INSTRUCTIONS FOR USE
3.		PROTECT FROM HEAT SOURCE AND RADIATION SOURCE
4.		KEEP DRY
5.		USE BY
6.		BATCH CODE
7.		CATALOGUE NUMBER
8.		STERILIZED USING ETHYLENE OXIDE
9.		DO NOT USE IF PACKAGE IS DAMAGED
10.		QUANTITY OF PRODUCT CONTAINED 1
11.		TEMPERATURE LIMITATION
12.		TYPE CF APPLIED PART
13.		DATE OF MANUFACTURE
14.		MANUFACTURER
15.		AUTHORISED REPRESENTATIVE IN THE EUROPEAN COMMUNITY

16.

**IMPORTER****【AFTER-SALES SERVICE】**

With “providing the medical sector with high quality and efficacy medical products” as its top operational objective, Shanghai MicroPort EP MedTech Co., Ltd. (hereinafter referred to as MicroPort EP) guarantees that its products are free of defects in materials or manufacturing when the clients receive them. For other questions relating to the products, please directly consult the company.

【SOLEMN STATEMENT】

Shanghai MicroPort EP MedTech Co., Ltd. (hereinafter referred to as MicroPort EP Co.) expressly states herein that its **EasyFinder™ Electrophysiology Diagnostic Catheter** is disposable and cannot be reused. MicroPort EP Co. will not recommend, indicate and imply in any manner the reusability of the system, and will not assume the responsibility for any accident or product damage resulting from reuse. **EasyFinder™ Electrophysiology Diagnostic Catheter** can be connected and used only with the compatible devices specified herein, and MicroPort EP Co. will not assume the responsibility for damage to product device, procedure failure and the like resulting from operating mistakes or any other human factors.

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Shanghai MicroPort EP MedTech Co., Ltd.

Address: Building 23&26&28, Lane 588, Tianxiong Rd., 201318 Shanghai, PEOPLE’S REPUBLIC OF CHINA

Postcode: 201318

Tel: +86 (21) 60969600

Fax: +86 (21) 20903925

E-mail: customerservice@everpace.com

Website: www.everpace.com

【AUTHORIZED REPRESENTATIVE IN THE EUROPEAN COMMUNITY】

MicroPort Medical B.V.

Address: Paasheuvelweg 25, 1105BP Amsterdam, The Netherlands

E-mail: cs@microport-int.com

Tel: +31(0)205450100-8

Fax: +31(0)205450109

**The copyright of this document is owned by Shanghai MicroPort EP MedTech Co., Ltd.
Copying is not allowed without its consent.**

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