Forlnk™

使用说明书	中文	页码	1
FORLNK™ 连接线缆			
Instructions for Use	English	Page	4
FORLNK™ Cable for CryoAblation Catheter			
Kullanım Talimatları	Türkçe	Sayfa	8
CryoAblasyon Kateteri için FORLNK™ Kablosu			
Οδηγίες χρήσης	Ελληνικά	Σελίδα	12
Καλώδιο FORLNK™ για καθετήρα CryoAblation			
Instrucciones de uso	Español	Página	16
Cable FORLNK™ para catéter de crioablación			
Instruções de utilização	Português	Página	20
Cabo FORLNK™ para cateter CrioAblação			
Ръководство за употреба	Български	страница	24
Кабел FORLNK™ за катетър за криоаблация			
Инструкции по применению	Русский	Страницы	28
Кабель FORLNK™ для катетера CryoAblation			
Bedienungsanleitung	Deutsch	Seite	32
FORLNK™ Kabel für CryoAblation-Katheter			
Instructions d'utilisation	Français	Page	36
Câble FORLNK™ pour cathéter de cryoablation			
Istruzioni per l'uso	Italiano	Pagina	40
Cavo FORLNK™ per catetere per crioablazione			
Gebruiksaanwijzing	Nederlands	Pagina	44
FORLNK™-kabel voor Cryo-ablatiekatheter			

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FORLNK™ Cable for CryoAblation Catheter

Shanghai MicroPort EP MedTech Co., Ltd.

【CABLE DESCRIPTION】

The cable is designed for the electric connection of IceMagic™ CryoAblation Catheter with IceMagic™ CryoAblation Console which are produced by Shanghai MicroPort EP MedTech Co., Ltd. (hereinafter called “Microport EP” for short). The cable may also be reused subjected to the reuse instructions provided below.

【PRODUCT SPECIFICATION】

Table 1: Specification of FORLNK™ Cable for CryoAblation Catheter

Cable type	Characteristic of Connector	Length
EPAU 250	34pins	2.5m

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The cable is designed for the electric connection of IceMagic™ CryoAblation Catheter with IceMagic™ CryoAblation Console.

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The cable should only be used by physicians and teams who are trained in electrophysiological techniques.

【WARNING】

- Use the cable prior to the “Use by” date on the package label.
- Patient or operator injury can result from improper handling of the cable.
- If the cable is used in the presence of electric equipment, noise may be induced into the cable.
- If conventional line powered equipment is used near the patient, follow electromechanical device and hospital standards to prevent extraneous electrical currents from reaching the electrophysiological equipment, catheter and heart. Such currents cause lethal arrhythmias.
- The cable should be visually inspect for any damage (i.e. kinks, nicks, crushed or elongated sections) prior to use or reuse. Do not attempt to repair any damages. In case of doubt, discard the cable and do not reuse.
- When used, cable shall be disposed of as medical waste according local laws and regulations.

【OPERATING INSTRUCTIONS】

Connection

- 1) Connect the cable to the handle of the IceMagic™ CryoAblation Catheter while maintaining sterility of the IceMagic™ CryoAblation Catheter.
- 2) Connect the cable to the IceMagic™ CryoAblation Console.

NOTE: The end of the connections should remain dry.

NOTE: Connecting cable to the IceMagic™ CryoAblation Console should be performed by personnel outside the sterile field.

Disconnection

- 1) Disconnect the catheter side of the cable: hold the catheter handle with one hand, while holding the arrow mark on the cable connector with the other hand and you can pull it out easily. Do not pull the body of the cable.
- 2) Disconnect the cryoablation console side of the: hold the arrow mark of the connector plug of the cable which is connected to the cryoablation console and pull it out.

【GUIDE FOR STERILIZATION AND REPROCESSING】

Please refer to the table 2 as below for the reprocessing instructions.

Table 2: Reprocessing Instructions of FORLNK™ Cable for CryoAblation Catheter

Manufacturer: Shanghai MicroPort EP MedTech Co., Ltd.	
Device: FORLNK™ Cable for CryoAblation Catheter	
WARNINGS	• Do not immerse the cable connectors in fluids. • Do not reuse the cable after ten (10) uses. • The cable should be visually inspected for any damage (i.e. cuts, kinks, nicks, crushed or elongated sections) prior to use or reuse. • Do not attempt to repair any damage. In case of doubt, discard the cable and do not use. • Do not expose cable to organic solvents. • Patient or operator injury can result from improper handling or storage of the cable. To prevent damage, store the cable in a clean and secure area
Limitations on processing	The cable is designed for a maximum of 10 processing cycles.
INSTRUCTIONS	
Initial treatment at the point of use	There are no requirements for initial treatment at point of use. The cable has been sterilized by Ethylene Oxide (ETO). The sterile state is effective only when the package is intact.
Preparation before cleaning	Visually inspect the cable for any damage such as crushed or elongated sections, cuts, kinks, or nicks. If damage is observed, do not use the cable and dispose the cable as medical waste according local laws and regulations .

Cleaning: Automated	Automated cleaning of the cable is not recommended because fluid may enter the connectors and leave moisture or residue in or on the connectors. Resulting adverse effects may include improper functioning (such as noisy signals or signal degradation) or arcing of electricity between the connector pins, which can lead to patient injury
Cleaning: Manual	- For each alligator clip: Pull off the protective sleeve in direction of the junction. - Clean and rinse with the following agents: - Cleaning agents with a pH-value 4....10 (note manufacturer's data) - Distilled water - Permissible cleaning method: - Manual cleaning by brushing down - Application time: 25 min - Cleaning Solution: RUHOF ENDOZIME® AW Plus W/APA - Concentration: 1:270 - Temperature: 40-45°C - Rinse with distilled water for 3min.
Disinfection: Automated	Automated disinfection of the cable is not recommended because fluid may enter the connectors and leave moisture or residue in or on the connectors. Resulting adverse effects may include improper functioning (such as noisy signals or signal degradation) or arcing of electricity between the connector pins, which can lead to patient injury
Disinfection: Manual	- Application time: note manufacturer's data, yet do not soak more than 1 h. Disinfect with the following agent: - Mixture of 70% isopropanol and 30% water - Immersion time: 25-30 min
Drying	Dry the cable carefully at a maximum of 120°C. Drying time: 25-30min
Maintenance, inspection and testing	Before packaging, visually inspect the cable. If visible soil is present, repeat the reprocessing instructions in this table. If there is evidence of damage, do not use the cable and dispose the cable as medical waste according local laws and regulations .
Packaging	Coil the cable into a loop when placing it into the pouch. Do not fold the cable because doing so may break the cable. Use a pouch large enough to prevent stress to the pouch seams and to prevent excessive bending of the cable.
Sterilization	- Pack the cable in the applicable sterile container, according to the respective sterilization procedure. - Sterilize with steam: 5min hold at 134°C
Storage	Follow permissible environmental conditions for storage.
Additional information	- The used product may be properly disposed of in contaminated hospital medical waste Never reuse the cable for more than 10 times. Reuse beyond this limit will compromise the normal function of the connection cable. MicroPort EP Co. will not assume the responsibility for the direct or subsequent adverse consequences resulting from reuse of the connection cable beyond the limit.
Manufacturer contact	For any questions relating to the cable, please directly consult the company. Please refer to the part of 【MANUFACTURER】 in this file for detailed contact information.

Note: The instructions provided above have been validated by the medical device manufacturer as being capable of preparing a medical device for reuse. It remains the responsibility of the processor to ensure that the processing as performed (using equipment, materials, and personnel in the processing facility) achieves the desired result. This requires validation and routine monitoring of the process. Likewise, any deviation by the processor from the instructions provided should be properly evaluated for effectiveness and potential adverse consequences.

【TRANSPORTATION REQUIREMENTS】

In transit, the product shall be protected from heavy load, direct sunlight and rain or as specified in the ordering contract.

【STORAGE REQUIREMENTS】

The product shall be stored in a shady and cool, dry, clean, and well-ventilated environment free of corrosive gases. The temperature during storage shall be kept in between 0°C and 45°C.

【SHELF LIFE】

While meeting the conditions for storage, the shelf life is two years.

Use this product before the end of the shelf life indicated on the package label.

【STERILIZATION】

This product has been sterilized with ethylene oxide gas. Do not use the product if the package is open or damaged. Use the product prior to the expiration date shown on the package label.

【SYMBOLS EXPLANATION】

1.		Consult Instructions for Use
2.		Protect from heat source and radiation source
3.		Keep dry
4.		Use by
5.		Batch Code
6.		Catalogue Number
7.		Sterilized using ethylene oxide
8.		Sterilized using steam or dry heat
9.		Do not use if package is damaged
10.		Quantity of product contained 1
11.		Temperature limitation
12.		Date of manufacture
13.		Manufacturer
14.		Authorized representative in the European Community
15.		Medical device
16.		Single sterile barrier system with protective packaging outside
17.		Importer
18.		UNIQUE DEVICE IDENTIFIER

【AFTER-SALES SERVICE】

With “providing the medical sector with high quality and efficacy medical products” as its top operational objective, Shanghai MicroPort EP MedTech Co., Ltd. (hereinafter referred to as MicroPort EP) guarantees that its products are free of defects in materials or manufacturing when the clients receive them. For other questions relating to the products, please directly consult the company.

【SOLEMN STATEMENT】

Shanghai MicroPort EP MedTech Co., Ltd. indicates definitely that the cable is used only for connecting IceMagic™ CryoAblation Catheter with IceMagic™ CryoAblation Console produced by Microport EP.

There is no proof to show that the accessory cable can connect to other brands of catheters and cryoablation consoles. Do not attempt to connect the accessory cable to other catheters or RF generator than IFU. Microport EP will not to be liable for any incidental, special, or consequential loss, damage, or expense resulting, directly or indirectly, from inaccurate connecting or any other man-caused accident.

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Shanghai MicroPort EP MedTech Co., Ltd.

Address: Building 23&26&28, Lane 588, Tianxiong Rd., 201318 Shanghai, PEOPLE'S REPUBLIC OF CHINA

Postcode: 201318

Tel: +86 (21) 60969600

Fax: +86 (21) 20903925

E-mail: customerservice@everpace.com

Website: www.everpace.com

【AUTHORIZED REPRESENTATIVE IN THE EUROPEAN COMMUNITY】

MicroPort Medical B.V.

Address: Paasheuvelweg 25, 1105BP Amsterdam, The Netherlands

E-mail: cs@microport-int.com

Tel: +31(0)205450100-8

Fax: +31(0)205450109

**The copyright of this document is owned by Shanghai MicroPort EP MedTech Co., Ltd.
Copying is not allowed without its consent.**

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**De copyright van dit document is eigendom van Shanghai MicroPort EP MedTech Co., Ltd.