使用说明书	中文	页码	1
标测导管连接尾线			
Instructions for Use	English	Page	4
FORLNK™ Cable for Mapping Catheter			
Kullanım Talimatları	Türkçe	Sayfa	8
Haritalama Kateteri için FORLNK™ Kablosu			
Οδηγίες χρήσης	Ελληνικά	Σελίδα	12
Καλώδιο FORLNK™ για Καθετήρα χαρτογράφησης			
Instrucciones de uso	Español	Página	17
Cable FORLNK™ para catéter de mapeo diagnóstico			
Instruções de Utilização	Português	Página	21
Cabo FORLNK™ para cateter de diagnóstico			
Инструкция за употреба	Български	Страница	25
Кабел FORLNK™ за катетър за картографиране			
Инструкция по применению	Русский	Страница	30
Кабель FORLNK™ для картирующего катетера			
Gebrauchsanweisung	Deutsch	Seite	34
FORLNK™ Kabel für Mappingkatheter			
Instructions d'utilisation	Français	Page	39
Câble FORLNK™ pour cathéter de cartographie HD			
Istruzioni per l'uso	Italiano	Pagina	43
Cavo FORLNK™ per catetere per mappaggio			
Instructies voor het gebruik	Nederlands	Pagina	47
FORLNK™-kabel voor Mappingkatheter			

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Shanghai MicroPort EP MedTech Co., Ltd.

【CABLE DESCRIPTION】



The cable is designed to connect a Electrophysiological Mapping Catheter which is produced by Shanghai MicroPort EP MedTech Co., Ltd. (hereinafter called “Microport EP” for short) to the appropriate recording equipment.

The cable may also be reused subjected to the reuse instructions provided below.

【PRODUCT SPECIFICATION】

There is 4 common-used specifications of cable for mapping catheter according to the pins (detailed in Table 1). Pins of cable (or sub-line seat) marked with numbers D, 2, 3, and etc, which indicated the pin corresponding to the electrode number of diagnostic catheter. The electrodes numbered from catheter remote end.

Table 1: Common-used specification of FORLNK™ Cable for Mapping Catheter

Terms	Characteristic of pins (connect to Standard recording)	Length	
		subject	Sub-line
COG010	10pins 	2m	18cm
COG004	4pins 	2m	18cm
COG014	14pins 	2m	18cm
COG026	26pins 	1.7m	20-25cm

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The cable is designed to connect a Electrophysiological Mapping Catheter to the appropriate recording equipment.

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This product shall be connected and used by doctors trained by the manufacturer's technicians or by the manufacturer's technicians themselves.

【WARNING】

- Prior to use, carefully check the product shelf life on the package and do not use product with damaged package and beyond the shelf life.
- Incorrect operation may cause injury to the patient and operator.
- Use of the connection cable with other electric equipment present may cause electrical interference.
- If motor and equipment with conventional circuit are present near the patient, comply with the standard for safety of medical electrical equipment (IEC 60601-1 3rd) and hospital standards to prevent current from approaching to electrophysiological equipment, catheter and the heart. This current may lead to fatal arrhythmia.
- Prior to use, carefully check if the connection cable is damaged in appearance (e.g. entanglement, notches, deformation or extension in connection cable). Do not try to repair the damaged cable. In case of any doubt, directly discard the connection cable.
- The rated accessory voltage of the catheter and cable is 450V and make sure that the output voltage used by the equipment is below this limit.
- When placing the catheter and cable, make sure that they are not in contact with the patient or other cables. When they are not in service for the time being, keep them in a place isolated from the patient.
- When used, cable shall be disposed of as medical waste according local laws and regulations.

【OPERATING INSTRUCTIONS】

Connection

- 1) Plug connector of the cable connect to the catheter, and the pins of the cable connect to the appropriate stimulation and/or recording equipment.
- 2) Connect the cable plug connector with the arrow mark to the catheter handle connector into the corresponding key bit slot: Hold the catheter handle, and push the cable connectors into the catheter handle connector. When you hear a sound "keta", that means the two connectors has locked to each other.
- 3) Insert the pins on the basis of pin numbers order into the corresponding jacks on the junction box of the standard EP recording equipment.

Disconnection

- 1) Disconnect the catheter side of the cable: Holding the catheter handle; another hand holding the arrow mark of the cable connector simultaneously and you can pull it out easily. Do not pull the body of the cable.
- 2) Disconnect the standard electrophysiological recording equipment side of the cable: Holding down the junction box; another hand holding the pins of the cable simultaneously and you can pull it out easily.

【PRODUCT TROUBLE-SHOOTING】

The following points can be used as instructions for Loop system of the diagnostic catheter / cable / electrophysiological recording system or Stimulation broken off.

- Cannot recognize the electrophysiological signals. If the recorder has a large number of electrophysiological signal interference or interruption, the reason is that the cable is not connected properly.

To solve the problem of the connector and the socket pins, please deal with it in the following order:

- 1) Please verify that the settings of multi-channel electrophysiological recorder or stimulator are match the channel of catheter and cable.
- 2) Please check and verify that all the jacks in the correct position when working. Confirm the electrophysiological recorder or stimulator connecting intact.
- 3) Pull out the plug or core needle, then reconnect it. Pull out pins and socket slightly, if the feedback signal is generated intermittently, there may be something wrong with the electrophysiological recorder or stimulator of connector / adapter.
- 4) If the signals are back to normal after replacing cable of Microport EP, there may be something wrong with the cable.

【GUIDE FOR STERILIZATION AND REPROCESSING】

Please refer to the table 2 as below for the reprocessing instructions.

Table 2: Reprocessing Instructions of FORLNK™ Cable for Mapping Catheter

Manufacturer: Shanghai MicroPort EP MedTech Co., Ltd.	
Device: FORLNK™ Cable for Mapping Catheter	
WARNINGS	● Do not immerse the cable connectors in fluids. ● Do not reuse the cable after ten (10) uses. ● The cable should be visually inspected for any damage (i.e. cuts, kinks, nicks, crushed or elongated sections) prior to use or reuse. ● Do not attempt to repair any damage. In case of doubt, discard the cable and do not use. ● Do not expose cable to organic solvents. ● Patient or operator injury can result from improper handling or storage of the cable. To prevent damage, store the cable in a clean and secure area
Limitations on processing	The cable is designed for a maximum of 10 processing cycles.

INSTRUCTIONS	
Initial treatment at the point of use	There are no requirements for initial treatment at point of use. The cable has been sterilized by Ethylene Oxide (ETO). The sterile state is effective only when the package is intact.
Preparation before cleaning	Visually inspect the cable for any damage such as crushed or elongated sections, cuts, kinks, or nicks. If damage is observed, do not use the cable and dispose the cable as medical waste according local laws and regulations .
Cleaning: Automated	Automated cleaning of the cable is not recommended because fluid may enter the connectors and leave moisture or residue in or on the connectors. Resultingadverse effects may include improper functioning (such as noisy signals or signal degradation) or arcing of electricity between the connector pins, which can lead topatient injury
Cleaning: Manual	● For each alligator clip: Pull off the protective sleeve in direction of the junction. ● Clean and rinse with the following agents: - Cleaning agents with a pH-value 4.....10 (note manufacturer's data) - Distilled water ● Permissible cleaning method: - Manual cleaning by brushing down - Application time: 25 min - Cleaning Solution: RUHOF ENDOZIME® AW Plus W/APA - Concentration: 1:270 - Temperature: 40-45°C ● Rinse with distilled water for 3min.
Disinfection: Automated	Automated disinfection of the cable is not recommended because fluid may enter the connectors and leave moisture or residue in or on the connectors. Resultingadverse effects may include improper functioning (such as noisy signals or signal degradation) or arcing of electricity between the connector pins, which can lead topatient injury

Disinfection: Manual	● Application time: note manufacturer's data, yet do not soak more than 1 h. Disinfect with the following agent: - Mixture of 70% isopropanol and 30% water - Immersion time: 25-30 min
Drying	Dry the cable carefully at a maximum of 120°C. Drying time: 25-30 min
Maintenance, inspection and testing	Before packaging, visually inspect the cable. If visible soil is present, repeat the reprocessing instructions in this table. If there is evidence of damage, do not use the cable and dispose the cable as medical waste according to local laws and regulations.
Packaging	Coil the cable into a loop when placing it into the pouch. Do not fold the cable because doing so may break the cable. Use a pouch large enough to prevent stress to the pouch seams and to prevent excessive bending of the cable.
Sterilization	● Pack the cable in the applicable sterile container, according to the respective sterilization procedure. ● Sterilize with steam: - 5 min hold at 134°C
Storage	Follow permissible environmental conditions for storage.
Additional information	● The used product may be properly disposed of in contaminated hospital medical waste ● Never reuse the cable for more than 10 times. Reuse beyond this limit will compromise the normal function of the connection cable. MicroPort EP Co. will not assume the responsibility for the direct or subsequent adverse consequences resulting from reuse of the connection cable beyond the limit.
Manufacturer contact	For any questions relating to the cable, please directly consult the company. Please refer to the part of 【MANUFACTURER】 in this file for detailed contact information.

Note: The instructions provided above have been validated by the medical device manufacturer as being capable of preparing a medical device for reuse. It remains the responsibility of the processor to ensure that the processing as performed (using equipment, materials, and personnel in the processing facility) achieves the desired result. This requires validation and routine monitoring of the process. Likewise, any deviation by the processor from the instructions provided should be properly evaluated for effectiveness and potential adverse consequences.

【TRANSPORTATION REQUIREMENTS】

In transit, the product shall be protected from heavy load, direct sunlight and rain or as specified in the ordering contract.

【STORAGE REQUIREMENTS】

The product shall be stored in a shady and cool, dry, clean, and well-ventilated environment free of corrosive gases. The temperature during storage shall be kept in between 0°C and 45°C.

【SHELF LIFE】

While meeting the conditions for storage, the shelf life is three years.

Use this product before the end of the shelf life indicated on the package label.

【STERILIZATION】

This product has been sterilized with ethylene oxide gas. Do not use the product if the package is open or damaged. Use the product prior to the expiration date shown on the package label.

【SYMBOLS EXPLANATION】

- | | | |
|----|---|---|
| 1. |  | Consult Instructions for Use |
| 2. |  | Protect from heat source and radiation source |
| 3. |  | Keep dry |
| 4. |  | Use by |
| 5. |  | Batch Code |
| 6. |  | Catalogue Number |
| 7. |  | Sterilized using ethylene oxide |
| 8. |  | Sterilized using steam or dry heat |

9.		Do not use if package is damaged
10.		Quantity of product contained 1
11.		Temperature limitation
12.		Date of manufacture
13.		Manufacturer
14.		Authorized representative in the European Community
15.		Medical device
16.		Single sterile barrier system with protective packaging outside
17.		IMPORTER
18.		UNIQUE DEVICE IDENTIFIER

【AFTER-SALES SERVICE】

With “providing the medical sector with high quality and efficacy medical products” as its top operational objective, Shanghai MicroPort EP MedTech Co., Ltd. (hereinafter referred to as MicroPort EP) guarantees that its products are free of defects in materials or manufacturing when the clients receive them. For other questions relating to the products, please directly consult the company.

【SOLEMN STATEMENT】

Shanghai MicroPort EP MedTech Co., Ltd. clearly declares: This product is used only for the electric connection of EasyFinder™ Fixed/Steerable curve diagnostic catheter or EasyLoop™ circular mapping catheter to the appropriate stimulation and/or recording equipment.

So far, there is no proof that this connection cable can be used with electrophysiological catheters or equipment of other brands. Do not use this product in the connection with catheters or RF generators beyond the scope specified by the operating instructions. Otherwise, MicroPort EP Co. will not assume the responsibility for the damaged product and equipment and surgery failures arising therefrom.

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Shanghai MicroPort EP MedTech Co., Ltd.

Address: Building 23&26&28, Lane 588, Tianxiong Rd., 201318 Shanghai, PEOPLE’S REPUBLIC OF CHINA

Postcode: 201318

Tel: +86 (21) 60969600

Fax: +86 (21) 20903925

E-mail: customerservice@everpace.com

Website: www.everpace.com

【AUTHORIZED REPRESENTATIVE IN THE EUROPEAN COMMUNITY】

MicroPort Medical B.V.

Address: Paasheuvelweg 25, 1105BP Amsterdam, The Netherlands

E-mail: cs@microport-int.com

Tel: +31(0)205450100-8

Fax: +31(0)205450109

**The copyright of this document is owned by Shanghai MicroPort EP MedTech Co., Ltd.
Copying is not allowed without its consent.**

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**De copyright van dit document is eigendom van Shanghai MicroPort EP MedTech Co., Ltd.